

PRODUKTINFORMATION

SAF-T HOLDER®



HÅLLARE TILL VENPROVTAGNINGSET

Vi kan erbjuda ett stort utbud av hållare i olika utföranden som används för venprovtagningar.

Hållare för venprovtagningsrör och olika modeller av blododlingsflaskor.



Benämning	Adapter	Specifikation	Art. nr	Enhet	Avdfp
Saf-T Holder®	Han luer	För vakumrör	96000	st	200
Saf-T Holder®	Hon luer	För vakumrör	96002	st	200
Saf-T Holder®	Han luer	För blododling	96004	st	50
Saf-T Holder®	Hon luer	För blododling	96005	st	50
Saf-T Holder®	*	Adapter pediatrika rör	96007	st	100
Saf-T Holder®	Han luer	För användning med Saf-T Wing®	96009	st	100

Smiths Medical Sverige AB

Box 1143, 164 22 Kista
Tlf. 08 - 594 772 50
Fax.08 - 594 772 69
info.sweden@smiths-medical.com
www.smiths-medical.com

Smiths Medical is part of the global technology business Smiths Group plc. Product(s)described may not be licensed or available for sale in all countries. The Smiths Medical design mark is a trademark of Smiths Medical. The symbol® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries.
©2012 Smiths Medical. All rights reserved.SE_11dGPI_STH-1