

HEINE BETA® 200 F.O. ENT Diagnostic Sets

 **NEW!** Now also available
with LED^{HD} illumination






Set complete with: **BETA200 F.O. Otoscope in XHL or LED**

1 set (4 pcs.) of reusable tips (B-000.11.111)

5 of both AllSpec disposable tips 2.5 and 4 mm Ø

spare bulb for XHL version

hard case

BETA 200 F.O. Otoscope Sets	2.5V XHL	3.5V XHL	LED
BETA battery handle	B-141.10.118		
 BETA4 USB rechargeable handle with USB cord and plug-in power supply		B-141.27.388	B-141.28.388
 BETA4 USB rechargeable handle		B-141.27.387	B-141.28.387
 BETA4 NT rechargeable handle with NT 4 table charger		B-141.23.420	B-141.24.420