

## HEINE mini 3000® Ophthalmoscope

### Compact Pocket Ophthalmoscope with XHL illumination



➤ Modern pocket instrument

➤ 5 different apertures

➤ Dust resistant

Modern pocket ophthalmoscope with XHL Xenon Halogen bulb technology. Concentrated bright light for perfect illumination, a brilliant image and accurate colour rendering. Fixation star aperture with grey filter for reducing the light intensity. Excellent optics and best colour rendering in minimum dimensions, ideal for rounds and house calls. Can be used with the mini3000 handle system. Available in black and blue.

- **HEINE optics produce crisp, clear images of the whole area under examination** and minimise reflections.
- **5 different apertures.** A full-function instrument incl. fixation star.
- Grey filter for light-sensitive patients.
- **Range of lenses: 18 lenses from -20 D to +20 D.**

+	1	2	3	4	6	8	10	15	20D
-	1	2	3	4	6	8	10	15	20D
- **Easy workflow**, because ophthalmoscopy is possible without dilation of the pupils.
- **Attachment clip with integrated on/off switch.** Switches off automatically when returned to the pocket.
- **Replaceable batteries.** Size AA, or optional rechargeable battery with mini NT table charger.\*
- **2-piece handle and head.** Easy to maintain, flexible.
- **Award-winning design.** Maximum quality with minimum dimensions.

mini 3000 Ophthalmoscope	2.5V XHL
with mini3000 battery handle, with batteries	<a href="#">D-001.71.120</a>
with mini3000 rechargeable handle (please order mini NT table charger separately)	<a href="#">D-001.71.200</a>
without handle	<a href="#">D-001.71.105</a>
XHL Xenon Halogen spare bulb	<a href="#">X-001.88.106</a>

#### mini 3000 apertures



Fixation star, large spot, small spot, hemispherical spot, red-free filter