


# TRU-CLOSE<sup>®</sup> Thoracic Vent

For treatment of pneumothorax


- For treatment of spontaneous, traumatic, or iatrogenic simple pneumothorax
- The self-contained vent allows complete ambulation during resolution
- All the components necessary for immediate insertion and evacuation are contained in the procedure tray
- The syringe-activated drain port allows aspiration of small amounts of fluid




The red signal diaphragm deflects upwards when the catheter tip enters the pleural space. It then continues to fluctuate with respiration as the intrapleural pressure varies from  $> +2\text{mm Hg}$  to  $< -2\text{mm Hg}$ .

When the signal diaphragm stays in the down position and the patient is asymptomatic the pneumothorax may be resolved.


- For treatment of spontaneous, traumatic, or iatrogenic simple pneumothorax
- The self-contained vent allows complete ambulation during resolution
- All the components necessary for immediate insertion and evacuation are contained in the procedure tray
- The syringe-activated drain port allows aspiration of small amounts of fluid




Over-the-wire insertion using optional Insertion Cannula (for TV11-13 & TV11-13E only.)



Quick and easy insertion using trocar and red signal diaphragm.



Double-valved aspiration cannula allows immediate manual air evacuation.



Allows easy attachment to a suitable suction system.



**Procedure Tray contains:**

- Thoracic vent
  - Trocar
  - Aspiration cannula
  - Suction tubing set
  - 60cc syringe
  - 3cc syringe
  - Safety needle, 25G x 5/8"
  - Safety needle, 22G x 1 1/2"
  - Scalpel
  - Filter Straw
  - Gauze sponges (2)
  - Fenestrated drape
  - CSR wrap
  - Ampule lidocaine 1%, 5 ml
  - Occlusion Plug
  - Povidone Iodine swabstick
- Note: Lidocaine and Povidone Iodine are only available in products sold in the U.S.A.

Local Representative:

**TRU-CLOSE Thoracic Vent Procedure Tray**

U.S.A. Codes	Int'l Codes	Cath. O.D. (Fr)	Cath. Length (cm)
TV11-10		11	10
TV11-13	TV11-13E	11	13
TV13-10	TV13-10E	13	10
TV13-13	TV13-13E	13	13

**Over-the-Wire Insertion Cannula for TV11-13(E)**

Code	Units/box
OTW11-13	5

**Suction Tubing Set**

Code	Units/box
TVST	5

**Occlusion Plug**

Code	Units/box
TVOCP	5

**Ordering Information: Call 1-800-538-7374**

For more information about UreSil products, please contact our sales department.

Made in the U.S.A. by UreSil, LLC.  
 5418 W. Touhy Ave., Skokie, IL 60077  
 Tel: (800) 538-7374 Tel: (847) 982-0200 Fax: (847) 982-0106  
 UreSil, Ureflex and TRU-CLOSE are registered trademarks of UreSil, LLC.  
 U.S. patent no. 4,664,660; 4,944,724; 4,966,197; 5,356,386; 6,770,063  
 SF 07-07 TCTV  
 EU: Tudor Health Care, Nottingham, NG1 1GF, Tel: 0044(0) 115 912 4277

