

# VICRYL\* Rapide (Polyglactin 910)

## SUTURE

### Description:

Braided Coated	For easy handling and secure knot tying. For smooth passage through tissue and easy knot tie down.
Synthetic Absorbable	For minimal tissue reaction. By hydrolysis. Essentially complete in 42 days.

### Tensile strength: Post implantation

	Approximate % original strength remaining
--	---

5 days	50%
10-14 days	0%

VICRYL\* Rapide suture typically falls off 7-10 days post-operative or can be wiped off subsequently with gauze. Normally the removal of the suture is not required.

**Color:** Undyed or violet.




**Range:** 8-0 to 1 (USP). Supplied as needled sutures and ligatures.





**Indications:** VICRYL\* Rapide is intended for use in soft tissue approximation where only short term wound support is required and where the rapid absorption of the suture would be beneficial. Typical areas of use include:


- Skin closure, cuticular or subcuticular
- Episiotomies
- Closure of oral mucosa
- Conjunctival sutures in ophthalmic surgery


Full details in the Instruction for Use included in every package.



<b>1/2 Circle Taper Point Needles</b>						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 RB-1 17 mm	5-0	○	70		36	V2130H
	4-0	○	70		36	V2140H
	3-0	○	70		36	V2150H
 SH-1 22 mm	4-0	○	70		36	V2180H
	3-0	○	70		36	V2190H
	2-0	○	70		36	V62H
 SH 26 mm	4-0	○	70		36	V4150H
	3-0	○	70		36	V4180H
	2-0	○	70		36	V4170H
	0	○	70		36	V66H




<b>1/2 Circle Taper Point Needles (Heavier Wire Diameter)</b>						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 CT-2 26 mm	2-0	○	70		36	V55H
	0	○	70		36	V56H
 CT-1 36 mm	2-0	○	90		36	V9450H
	0	○	90		36	V9480H
	1	○	90		36	V70H
 CT 40 mm	0	○	90		36	V9580H
 CTX 48 mm	0	○	90		36	PN1695H

<b>1/2 Circle ETHIGUARD* Blunt Point Needles</b>						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 CTB-1 36 mm	2-0	○	90		36	W9952H


<b>Straight Taper Point Needles</b>						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 ST-4 19 mm	4-0	○	70		36	BA1673H

**VICRYL\* Rapide**  
(Polyglactin 910)


**1/2 Circle TAPERCUT\* Needles**

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 V-5 17 mm	3-0	○	70		12	OA1674G
 V-7 26 mm	3-0	○	70		12	OB1674G
	3-0	○	70		36	OB1674H
 V-34 36 mm	2-0	○	90		36	V5170H
	0	○	90		36	V5180H





**Compound Curve TAPERCUT\* Needles**

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 DA 25 mm	4-0	WHITE	45		12	V2912G






**1/2 Circle Reverse Cutting PRIME\* Needles**

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 P-2 8 mm	4-0	○	45		36	JA1642H




**1/2 Circle Reverse Cutting Needles**

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 J-1 17 mm	4-0	○	45		12	18510G
	4-0	○	45		36	V4550H
	3-0	○	45		36	V76H
 X-1 22 mm	3-0	○	45		36	V7902H
	2-0	○	45		36	V81H
 CP-1 36 mm	2-0	○	90		36	V4721H
 CP 40 mm	0	○	90		36	V4851H
	1	○	90		36	V4861H



**3/8 Circle Reverse Cutting PRIME\* Needles**


Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 P-1 11 mm	6-0	○	45	MULTIPASS*	36	MPVR9913H
	6-0	●	45		36	V32H
	6-0	○	45		12	W9913
	5-0	○	45	MULTIPASS*	36	MPVR9915H
	5-0	○	45		12	W9915
 P-3 13 mm	5-0	○	45	MULTIPASS*	36	MPVR4930H
	5-0	○	45		36	V4930H
	4-0	○	45	MULTIPASS*	36	MPVR4940H
	4-0	○	45		36	V4940H
 PS-2 19 mm	4-0	○	45	MULTIPASS*	36	MPVR4960H
	4-0	○	75		12	W9922
	4-0	○	75	MULTIPASS*	36	MPVR4260H
	3-0	○	45		36	V4970H
	3-0	○	45	MULTIPASS*	36	MPVR4970H
	3-0	○	70	MULTIPASS*	36	MPVR4270H
	3-0	○	75		12	W9923
 PS-2 S 19 mm	4-0	○	45	Slimmed	36	V4960H
	4-0	○	70	Slimmed	36	V4260H
 PS-1 24 mm	3-0	○	45		36	V6830H
	3-0	○	70	MULTIPASS*	36	MPVR9360H
	3-0	○	70		36	V9360H
	2-0	○	70	MULTIPASS*	36	MPVR9380H
	2-0	○	70		36	V9380H





**3/8 Circle Reverse Cutting Needles**

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 FS-2 19 mm	5-0	○	45		36	V2910H
	4-0	○	45		36	FH1642H
	3-0	○	45		12	V2930G
	3-0	○	45		36	V2930H
 FS-2S 19 mm	4-0	○	45	Slimmed	36	V2920H
 FS-1 24 mm	3-0	○	45		36	V4421H
	3-0	○	70		36	V4420H
	2-0	○	45		36	V4431H

**VICRYL\* Rapide**  
(Polyglactin 910)

▲ Straight Cutting Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 TS 48 mm	3-0	○	75		36	VR2239
	2-0	○	75		36	VR2240
 KS 60 mm	3-0	○	70		36	AA1673H

▲ 3/8 Circle Conventional Cutting PRIME* Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 PC-3 16 mm	4-0	○	70	MULTIPASS*	36	MPVR9918H
	3-0	○	75		12	W9919

▼ 3/8 Circle MICROPPOINT* Spatula						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 GS-12 5,5 mm	8-0	●	30		12	V1220G
 GS-9 6,5 mm	8-0	●	30		12	V1210G
	7-0	●	30		12	V1211G
 GS-9, GS-9 6,5 mm	8-0	●	30		12	V39G
	7-0	●	45		12	V38G
 TG140-8 6,5 mm	7-0	●	30		12	W9914

Ligatures						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
SUTUPAK*	4-0	○	3x45		24	V8643E
	2-0	○	3x45		24	V8645E
	0	○	3x45		24	V8646E
	0	○	6x45		24	V8980E