

MONOCRYL*

(Poliglecaprone 25) Monofilament

SUTURE

Description: **Monofilament** For smooth passage through tissue.
Synthetic For minimal tissue reaction.
Absorbable By hydrolysis. Essentially complete in 90 -120 days.

Tensile strength: **Post implantation** **Approximate % original**
7 days **60% (violet)**
50% (undyed)

14 days **30% (violet)**
20% (undyed)

21 days **0% (undyed)**
28 days **0% (violet)**

Color: Violet or undyed.



















Range: 6-0 to 1 (USP). Supplied as needed sutures and ligatures.

Indications: MONOCRYL* sutures are intended for use in general soft tissue approximation and/or ligation where an absorbable material is indicated. Typical areas of use include:





- Subcuticular sutures
- Small intestine anastomoses
- Urological anastomoses




Full details in the Instruction for Use included in every package.






 1/2 Circle Taper Point Needles							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 TF 13 mm	6-0	●	70		36	Y432H	
	5-0	●	70		36	Y925H	
 RB-1 17 mm	6-0	●	70		36	Y212H	
	5-0	●	70		36	Y213H	
	4-0	●	70		36	Y214H	
	3-0	○	70		36	T4215AH	
	3-0	●	70		36	Y215H	
 JRB-1  VISI-BLACK 17 mm	4-0	●	70	VISI-BLACK*	36	Y3040H	
 JRB-1 JRB-1  VISI-BLACK 17 mm	5-0	●	70	VISI-BLACK*	36	Y3213H	
	4-0	●	70	VISI-BLACK*	36	Y3214H	
 SH-1 22 mm	4-0	●	70		36	Y218H	
	3-0	●	70		36	Y219H	
 JB-1  VISI-BLACK 22 mm	4-0	●	70	VISI-BLACK*	36	Y3100H	
	3-0	●	70	VISI-BLACK*	36	Y3110H	
	2-0	●	70	VISI-BLACK*	36	Y3120H	
 JB-1 JB-1  VISI-BLACK 22 mm	4-0	●	70	VISI-BLACK*	36	Y3218H	
	4-0	●	90	VISI-BLACK*	36	Y5218H	
	3-0	●	70	VISI-BLACK*	36	Y3219H	
 SH 26 mm	4-0	●	70		36	Y415H	
	3-0	○	70		12	NP0001	
	3-0	●	70		36	Y416H	
	2-0	●	70		36	Y417H	
 JB  VISI-BLACK 26 mm	4-0	●	70	VISI-BLACK*	36	Y3150H	
	3-0	●	70	VISI-BLACK*	36	Y3160H	
	3-0	●	8x45	CONTROL RELEASE*, VISI-BLACK*	12	Y3864G	
	2-0	●	70	VISI-BLACK*	36	Y3170H	
 JB JB  VISI-BLACK 26 mm	3-0	●	70	VISI-BLACK*	36	Y3416H	
	3-0	●	90	VISI-BLACK*	36	Y5270H	
	2-0	●	70	VISI-BLACK*	36	Y3417H	
 MH 36 mm	2-0	●	70		36	Y227H	

MONOCRYL*
(Poliglecaprone 25) Monofilament







1/2 Circle Taper Point Needles (Heavier Wire Diameter)						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 CT-2 26 mm	2-0	●	90		36	Y326H
	0	●	70		36	Y270H
 CT-1 36 mm	2-0	●	70		36	Y250H
	0	●	70		36	Y260H
	0	●	90		36	Y946H
	1	●	70		36	Y261H
 CT 40 mm	1	●	90		36	Y950H
 CTX 48 mm	2-0	●	90		12	Y979G
	1	●	90		12	Y977G

5/8 Circle Taper Point Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 UR-6 26 mm	3-0	●	70		12	W3625
	2-0	●	70		36	Y605H
 UR-6 26 mm	3-0	●	70		36	Y602H
	2-0	●	70		36	A1H
 UR-5 36 mm	2-0	●	70		36	Y375H


1/2 Circle TAPERCUT* Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 V-7 26 mm	2-0	●	70		36	Y999H
 V-34 36 mm	2-0	●	70		12	W3488
	0	●	70		12	W3489
	1	●	90		36	Y919H

Compound Curve TAPERCUT* Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 DA-1 13 mm	5-0	○	45		12	Y4950G




3/8 Circle Reverse Cutting PRIME* Needles



Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 P-1	11 mm	6-0	45	MULTIPASS*	36	MPY489H
		6-0	45		36	Y489H
		5-0	45	MULTIPASS*	36	MPY490H
 P-3	13 mm	6-0	45	MULTIPASS*	36	MPY492H
		6-0	45		36	Y492H
		5-0	45	MULTIPASS*	36	MPY493H
		5-0	45		36	Y491H
		5-0	45		36	Y493H
		4-0	45	MULTIPASS*	36	MPY494H
		4-0	45		36	Y494H
 PS-3	16 mm	5-0	45	MULTIPASS*	36	MPY500H
		5-0	45		12	W3204
		4-0	45	MULTIPASS*	36	MPY501H
		4-0	45		12	W3205
 PS-2	19 mm	5-0	45	MULTIPASS*	36	MPY495H
		4-0	45	MULTIPASS*	36	MPY496H
		4-0	45		12	W3206
		4-0	70	MULTIPASS*	36	MPY426H
		3-0	45	MULTIPASS*	36	MPY497H
		3-0	45		12	W3207
		3-0	70		36	C427
		3-0	70	MULTIPASS*	36	MPY427H
 PS-1	24 mm	4-0	45		36	T4682H
		4-0	70	MULTIPASS*	36	MPY935H
		4-0	70		36	Y495H
		3-0	70		12	W3208
		3-0	70	MULTIPASS*	36	MPY683H
		2-0	70	MULTIPASS*	36	MPY684H
		2-0	70		36	Y684H
 PS	26 mm	4-0	70	MULTIPASS*	36	MPY3212H
		4-0	70		12	W3212
		3-0	70	MULTIPASS*	36	MPY3213H
		3-0	70		12	W3213


Compound Curve Reverse Cutting PRIME* Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 PS-4C	16 mm	5-0	45		12	NP0127

MONOCRYL*
(Poliglecaprone 25) Monofilament

▼ 3/8 Circle Reverse Cutting Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 FS-2 19 mm	4-0	○	45		36	Y292H
	4-0	○	70		36	Y422H
	3-0	○	45		36	Y293H
	3-0	○	70		36	Y423H
 FS-1 24 mm	4-0	○	70		36	Y291H
	3-0	○	70		36	Y442H
	2-0	○	70		36	Y443H
 FS 28 mm	3-0	○	70		12	W3326

▲ Straight Cutting Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 TS 48 mm	4-0	○	70		36	C1041
 KS 60 mm	4-0	○	70		12	W3649
	3-0	○	70		12	W3650

▼ 1/4 Circle Spatula Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 S-14 8 mm	6-0	○	45		12	W3552

Ligatures						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
LIGAPAK*	3-0	●	250	Reel	10	W3020
	2-0	●	250	Reel	10	W3021
SUTUPAK*	2-0	●	150		12	W3113
	1	●	150		12	W3115
	0	●	150		12	W3114
	2-0	●	2x70		36	Y625H
	3-0	○	3x45		36	Y634H