

**ETHILON\***  
Monofilament Polyamide 6

**SUTURE**  
Description:

<b>Monofilament</b>	For smooth passage through tissue.
<b>Synthetic</b>	For minimal tissue reaction.
<b>Non-absorbable</b>	Provides prolonged tensile strength retention in tissue.

**Color:** Blue, black or clear.








**Range:** 11-0 to 2 (USP). Supplied as needed sutures armed with, for example, MULTIPASS\* and PRIME\* needles. ETHILON\* is also available as ligatures.

**Indications:** ETHILON\* (black and clear) sutures are intended for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, microsurgical and neurosurgical procedures. Blue ETHILON\* sutures are intended for use in skin closure. Typical areas of use include:




- Cuticular sutures
- Nerve adaptation
- Ophthalmology









Full details in the Instruction for Use included in every package.












3/8 Circle Reverse Cutting PRIME* Needles							
Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 P-6	8 mm	7-0	●	45		36	1647H
		7-0	●	45		36	1696H
 P-1	11 mm	7-0	●	45	MULTIPASS*	36	MPE1696H
		6-0	●	45		36	697H
		6-0	●	45	MULTIPASS*	36	MPE697H
		6-0	●	45		36	1698H
 P-3	13 mm	6-0	●	45		36	1698H
		6-0	●	45	MULTIPASS*	36	MPE1698H
		5-0	●	45		36	698H
		5-0	●	45		36	EH7800H
		5-0	●	45	MULTIPASS*	36	MPE698H
		5-0	●	45	MULTIPASS*	36	MPE7800H
		5-0	●	45		24	W1852T
		4-0	○	45		36	691H
		4-0	●	45		36	699H
		4-0	●	45		36	EH7801H
		4-0	○	45	MULTIPASS*	36	MPE691H
		4-0	●	45	MULTIPASS*	36	MPE699H
		4-0	●	45	MULTIPASS*	36	MPE7801H
		 PS-3	16 mm	6-0	●	45	MULTIPASS*
6-0	●			45		24	W1614T
5-0	●			45		36	EH7930H
5-0	●			45	MULTIPASS*	36	MPE7930H
4-0	●			45		36	EH7931H
4-0	●			45	MULTIPASS*	36	MPE7931H
4-0	●			45		24	W1861T
 PS-2	19 mm	5-0	●	45		36	EH7804H
		5-0	●	45	MULTIPASS*	36	MPE7804H
		4-0	●	45		36	EH7162H
		4-0	●	45		24	W1620T
		4-0	●	45	MULTIPASS*	36	MPE1667H
		4-0	●	45	MULTIPASS*	36	MPE1855H
		3-0	●	45		36	EH7803H
		3-0	●	45	MULTIPASS*	36	MPE1621H
 PS-2 S	19 mm	4-0	●	45	SLIMMED	36	1667H
 PS-1	24 mm	4-0	●	45		36	EH7957H
		4-0	●	75		36	EH7958H
		4-0	●	75	MULTIPASS*	36	MPE7958H
		3-0	●	45		24	W1857T
		3-0	●	75		36	EH7954H
		3-0	●	75	MULTIPASS*	36	MPE7954H

**ETHILON\***  
Monofilament Polyamide 6





3/8 Circle Reverse Cutting PRIME* Needles							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 PS	26 mm	3-0	●	45	MULTIPASS*	36	MPE1625H
		3-0	●	45		24	W1625T
		3-0	●	75		24	W1685T
 PSL	30 mm	3-0	●	75		36	EH7113H
 PSLX	36 mm	3-0	●	75		24	W1627T


3/8 Circle Reverse Cutting Needles							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 G-2	13 mm	6-0	●	45		36	667H
		5-0	●	45		36	EH7945H
 M-2	13 mm	5-0	●	45		36	644H
 FS-3	16 mm	6-0	●	35		12	W506
		6-0	●	45	VISI-BLACK*	36	EH7177H
		5-0	●	45		36	1628H
		5-0	●	45		36	EH7823H
		4-0	●	45		36	EH7824H
 FS-2	19 mm	5-0	●	45		36	EH7790H
		4-0	●	45		36	EH7144H
		4-0	●	75		36	EH7143H
		3-0	●	45		36	EH7792H
 FS-2S	19 mm	4-0	○	45	SLIMMED	36	611H
		4-0	●	45	SLIMMED	36	EH7791H
 FS-1	24 mm	5-0	●	45		36	EH7640H
		4-0	●	45		36	EH7825H
		4-0	●	45		36	S2450
		4-0	●	75		36	FJ5272H
		3-0	●	45		36	EH7794H
		3-0	●	75		36	EH7795H
		2-0	●	45		36	EH7826H
		2-0	●	45		36	EH7826H
 FS	26 mm	4-0	●	45		36	FK5242H
		3-0	●	45		36	ES163H
		2-0	●	45		36	EH7796H
 FSL	30 mm	4-0	●	75		36	1670H
		3-0	●	45		36	EH7936H
		3-0	●	75		36	EH7936H
		2-0	●	75		36	EH7798H


3/8 Circle Reverse Cutting Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 FSLX 36 mm	3-0	●	75		36	1673H
	2-0	●	75		36	1674H
	0	●	75		36	EH7809H
 LR-40 60 mm	2-0	●	100		12	W558
 LR 76 mm	2	●	3x50		12	3460TG



3/8 Circle Conventional Cutting PRIME <sup>®</sup> Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 PC-1 13 mm	6-0	●	45		36	JM5140H
 PC-3 16 mm	6-0	●	45	MULTIPASS <sup>®</sup>	36	MPE1853H
	6-0	●	45		24	W1815T
	6-0	●	45		24	W1853T
	5-0	●	45	MULTIPASS <sup>®</sup>	36	MPE1816H
	5-0	●	45		24	W1816T
	5-0	●	45		24	W1854T
 PC-5 19 mm	4-0	●	45	MULTIPASS <sup>®</sup>	36	MPE1856H
	4-0	●	45		24	W1856T
 PC-25 26 mm	6-0	●	45		24	W1859T
	5-0	●	45		24	W1858T
	5-0	●	45		12	W532
 KP-3 90 mm	2-0	●	100		12	W795
	1	●	100		12	W797
	1	●	2x50	With retention tubing	12	W2797
 KP-3 90 mm	1	●	75		12	W689


**ETHILON\***  
Monofilament Polyamide 6

▼ 1/2 Circle Reverse Cutting Needles (Heavier Wire Diameter)						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 M-1	13 mm	5-0	●	45		36 648H
 J-1	17 mm	5-0	●	75		36 F2423
 X-1	22 mm	3-0	●	45		36 EH7943H
 CPX	48 mm	1	●	100		12 W738
















▲ Straight Cutting Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 KS	60 mm	3-0	●	75		36 627H
		2-0	●	75		36 628H

◎ 1/2 Circle Taper Point Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 MH-1	31 mm	0	●	150	Loop	12 W739




◎ 1/2 Circle Taper Point Needles (Heavier Wire Diameter)						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 CT	40 mm	0	●	200	Loop	12 W760
 CTX	48 mm	1	●	150	Loop	12 W748

◎ Straight Taper Point Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 STAW	65 mm	2-0	●	75		24 EH7625E




 3/8 Circle Taper Point Needles


Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
BV 140-3 	3,7 mm	10-0	●	15	Non-sterile	12 NS2820
BV130-3 	3,7 mm	9-0	●	13		12 FG2819
BV75-3 	3,8 mm	10-0	●	13		12 W2870
		9-0	●	13		12 W2871
BV50-3 	3,8 mm	11-0	●	13		12 W2881
BV130-4 	4,7 mm	9-0	●	13		12 W2813
		8-0	●	13		12 W2812
BVH100-3 	5 mm	10-0	●	13		12 W2880
BV75-4 	5,1 mm	10-0	●	13		12 W2850
BV-6 	5,1 mm	11-0	●	10		12 EH7438G
BV-4 	5,1 mm	10-0	●	15		12 2814G
		10-0	●	13		12 W2830
		9-0	●	13		12 W2829
BV100-4 	5,1 mm	9-0	●	15		12 EH7448G
		10-0	●	15		12 2810G
		9-0	●	15		12 2809G
BV-2 	6,5 mm	8-0	●	15		12 2808G
		10-0	●	13		12 W2810
		8-0	●	13		12 W2808
BV130-5 	6,5 mm	10-0	●	30		12 W2801
BV175-6 	8 mm	8-0	●	45		12 W8170
BV-1 	9,3 mm	7-0	●	45		36 EH7448H
C-1 	13 mm	6-0	●	75		36 EH7500H



 3/8 Circle CS ULTIMA\* Spatula


Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
UCS-35 	6,2 mm	10-0	●	15		12 7060G
CS-140-6 	6,5 mm	10-0	●	30		12 U7003
CS-140-6 	6,5 mm	10-0	●	30		12 U7061


**ETHILON\***  
Monofilament Polyamide 6


▼ 3/8 Circle MICROPOINT* Spatula						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
CS-35-C 	6 mm	10-0	●	15	12	EH7980G
GS-9 	6,5 mm	10-0	●	30	12	1718G
		9-0	●	30	12	1717G
TG140-6  TG140-6	6,5 mm	11-0	●	30	12	W1780
		10-0	●	30	12	W1770
		9-0	●	30	12	W1769



▼ 1/2 Circle MICROPOINT* Reverse Cutting Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
G-7 	8 mm	8-0	●	45	12	W1765


▼ 1/2 Circle MICROPOINT* Spatula						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
TG160-6 	5,5 mm	10-0	●	30	12	W1758
TG175-8 	7 mm	8-0	●	30	12	1716G

▼ 1/2 Circle CS ULTIMA* Spatula						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
CS-160-6 	5,5 mm	10-0	●	30	12	U7000

▼ 1/4 Circle MICROPOINT* Spatula						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
TG100-8 	6,5 mm	8-0	●	30	12	B7714

▼ 1/4 Circle Spatula Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
S-14 	8 mm	5-0	●	20	36	MN5132H

1/4 Circle TAPERCUT* Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 V100-3      3,5 mm	10-0	●	13		12	W2899
	9-0	●	13		12	W2898
 V150-4      4,5 mm	10-0	●	15		12	W1701

Compound Curve CS ULTIMA* Spatula						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 CS-C-6      5 mm	10-0	●	30		12	U7009

Ligatures						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
SUTUPAK*	4-0	●	3x45		36	EH6751H
	3-0	●	10x45		36	EH6562H
LIGAPAK*	2-0	●	250	Reel	20	F4354