

TravelJohn Urine collector



This unique urine collector, that consists of a resealable outer pouch, and an inner pouch with a superabsorbent powder, that changes the urine to a non smelling gel.

TravelJohn is used instead of the wellknown plastique urine collectors, that often brings the problems with spill off, wich means that the urine will be either in the bed or on the floor. A normal plastique urine collector often has a unpleasant smell.

How does TravelJohn work

The superabsorbent powder in the inner pouch absorbs the urine and changes that to a gel. This very special Liqisorb gel has an antimicrobial feature, that neutralizes bacteria and smell.

The pouch is biodegradable, that mean the pouch can be destroyed as normal waste.

With TravelJohn Urine collector you secure

- Urine is safely collected in the bag
- Urine is changed to a consistant gel
- The absorbent powder destroys smell and bacteria
- The pouch can hold up to 800 ml
- The pouch is resealable, which means that it can be used more than one time

Where do you use the TravelJohn Urine Collector -Health care

- In hospitals, secures no spillage and smell
- In nursing homes
- Intermittent chatheter users with no access to a toilet
- People in wheelchairs, primarily men

Private users

- Saylor
- Camping
- Festivals
- Traveling