

PRODUCT SPECIFICATION

REF 1925 BREATHING MASK IN KEYRING POUCH



1 GENERAL PRODUCT DESCRIPTION

Description

A breathing mask that consists of a one-way valve (green upper part, red lower part and a membrane) and a plastic sheet with graphic instructions printed on it. The Breathing mask is packed in a small pouch, which can be used as a key ring. The pouch is delivered in a plastic zip bag.



The device is intended to be used as an aid of assistance when giving mouth-to-mouth resuscitation to an injured person. The device is also intended to work as a barrier and limit the contact between the resuscitator and the injured person's blood, secrete, vomit and other body fluids.

CE-labelled.

Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.

2 DETAILED PRODUCT INFORMATION

Content per consumer unit

<u>Product/Material</u>	<u>Quantity</u>	<u>Dimensions (mm)</u>
Breathing mask (plastic sheet with valve)	1	27 ± 0,5 x 25 ± 0,5 Thickness 0,1 mm

Material

Plastic sheet	PVC with phthalate (DEHP, DMP, BBP, DINP, DIDP, DNOP) contents < 0.003%
Valve, green upper part	Green ABS- plastic
Valve, red lower part	Red ABS-plastic
Valve, sealing	Red ABS-plastic
Valve, membrane	Red TPE (Thermoplastic polyester)
Plastic zip bag	Low density polyethylene

Shelf life: 5 years from production date

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310611925016	7310801925000	7312020012196
Content:	1 piece	40 CU	15 SKU/layer / 5 layers
Size (w x l x h mm):	80 x 44 x 137	205 x 260 x 180	Pallet 800 x 1200x 1050
Gross weight:	20 g	976 g	98 kg
Packing material:	Plastic 2 g	Carton 155 g	-

