

Vario 18



VERSATILITY AND EFFICIENCY



Vario 18 Suction Pump

Vario 18

Versatile and efficient suction

Vario 18 is an efficient medical suction pump. An all-rounder for hospitals, clinics and homecare, it offers reliability, low noise and mobility. Its innovative QuatroFlex drive unit allows for easy servicing.

Vario 18 can be used with either a Disposable or a Reusable Collection System, according to each customer's needs and demands.

Manufactured with swiss precision and quality, **Vario 18** offers a range of useful features:



- a Overflow protection/bacteria filter
- b Vacuum gauge in kPa and mmHg
- c LED warning light
- d Carrying handle
- e Membrane vacuum regulator
- f Tubing holder
- g Plug for 12 VDC IN car connection cable

Unique swiss technology

The innovative 4-fold piston/cylinder system, QuatroFlex, builds up the vacuum within seconds. It also convinces through compact size and low noise during operation. As well as being unique in design and technology it is easy to service and maintain. A fact, highly appreciated by medical technicians. The membrane vacuum regulator allows for precise vacuum settings. It needs to be pressed down to change the vacuum, thus making it impossible to change the settings inadvertently. The Medela Vario 18 can therefore also be used at the patient's bedside.



Precise vacuum setting with press-down knob

Versatility for a wide range of applications

Vario 18 is lightweight and portable for use across departments in hospitals and clinics and at home. It can be used for small surgical interventions as well as for airway suctioning. Variable suction levels can conveniently be set with the vacuum regulator. The Vario 18 can be equipped with either a Disposable or a Reuseable Collection System.



Optimal hygiene for many applications

Hygienic and easy use

Vario 18 offers easy, intuitive handling. All accessories are designed to complement it perfectly and to prevent false connections. Optional accessories include a trolley, carrying bag and various attachments. It also features an inbuilt battery charger that allows for mobile use. The need to recharge is made audible and visible by special battery alarms.



Can be equipped for individual needs

Fast Facts

- Medela offers healthcare professionals 50 years of expertise in medical vacuum technology.
- Our swiss development and manufacturing process yields mechanical precision to guarantee excellent performance.
- QuatroFlex technology: silent function: can be used in environments in which quiet operation is necessary.
- Precise vacuum settings and fast vacuum build-up: ready when needed.
- One device for multiple needs.
- New users need little instruction to operate.
- Mains-independent for 30 minutes.
- Acoustic and optical battery management.
- Wide range of accessories available for individual application.
- Vario 18 incorporates reliability and ease of use like the Medela suction pumps for operating theatre.
- 2 years warranty on device, 6 months warranty on battery/QuatroFlex.

Accessories

Filters

to protect the pump from overflow, the environment from bacteria and to neutralise odour.



Trolley

with rails for mobile transport of the pump and accessories.



Reusable Collection System

autoclavable jars and lids collect secretions easily and economically.
 Sizes available:
 1, 2, 3 and 5 litres.



Holders and plates

can be used to attach Vario 18 to standard rails.



Disposable Collection System

with disposable suction liners and reusable PC jars that offer a hygienic, simple handling and efficient fluid management.
 Sizes available:
 1.5 and 2.5 litres.



Carrying bag

for transporting the pump.





Use Medela accessories for even greater flexibility and ease of use.



Plug for 12 VDC IN car connection cable




for operating the pump in the car.



Technical data

	High vacuum -75 kPa -563 mmHg
	low flow 18 L/min.

	AC 3,5 kg AC/DC 4,2 kg
	H x W x D 380 x 170 x 285 mm

	IPX1
	ISO 9001 ISO 13485 CE (93/42/EEC), IIa 2 years warranty
	

Medical Vacuum Technology for Healthcare Professionals

Please contact us or your local Medela representative for details.

 Medela AG
 Lättichstrasse 4b
 6341 Baar / Switzerland
www.medela.com

Local contact:

Medela Australia Pty Ltd
 Medical Technology
 3 Arco Lane
 Heatherton, Vic, 3202
 Australia
 Phone +61 (0) 3 9552 8600
 Fax +61 (0) 3 9552 8699
contact@medela.com.au
www.medela.com.au

Medela Canada Inc.
 4160 Sladeview Cres., Unit #8
 Mississauga, Ontario
 L5L 0A1
 Canada
 Toll Free 800.435.8316
 Toll Free Fax 800.995.7867
info@medela.ca
www.medela.ca

Medela UK
 Huntsman Drive
 Northbank Industrial Park
 Irlam, Manchester M44 5EG
 UK
 Phone +44 870 950 5994
 Fax +44 870 389 2233
info@medela.co.uk
www.medela.co.uk

Medela Inc. USA
 P.O. Box 660
 1101 Corporate Drive
 McHenry, IL 60050
 USA
 Phone +1 877 735 1626
 Fax +1 815 363 2487
suction@medela.com
www.medelasuction.com