


Tubings


For hygienic transportation of secretions, reliable tubings complete the high quality system. Depending on the application, Medela offers various materials, diameters and connection pieces with the tubing.

Silicone tubings



Ref. 077.0901 (1 m)
 Ref. 077.0902 (2 m)
 Ref. 077.0905 (5 m)
 Ref. 077.0900 (25 m)

Silicone tubings
 – Ø 7 x 12 mm
 – Lengths 1, 2, 5, or 25 m



Ref. 077.0912 (60 cm)
 Ref. 077.0913 (100 cm)

Silicone tubings with 1 coupling piece
 – Ø 7 x 12 mm
 – Lengths 60 or 100 cm



Ref. 077.0922 (60 cm)
 Ref. 077.0921 (30 cm)

Silicone tubings with 2 coupling pieces
 – Ø 7 x 12 mm
 – Lengths 30 or 60 cm



Ref. 077.0931
 Ref. 077.0185 (60 cm)

Silicone tubing with 1 coupling piece and 1 coupling piece 90° angle
 – Ø 7 x 12 mm
 – Length 25cm



Ref. 077.0970

Silicone tubing
 – Ø 5 x 10 mm
 – Length 25 m

Disposable tubing



Ref. 077.0951

Disposable tubing with fingertip
 – PVC, sterile
 – Ø 6 x 9 mm
 – Length 180 cm



Ref. 077.0965
 Ref. 077.0952
 Ref. 077.0191

Disposable tubing with funnel connectors, Ø 5.5 x 8.5 mm
 – with two funnel connections, 150 cm
 – with one funnel connection (non-sterile), 180 cm
 – with one funnel connection, 300 cm



Ref. 077.0184

Disposable tubing with universal connectors
 – for use with Single Use Bird Cup or reusable Bird Cups
 – PVC, sterile
 – Ø 5.5 x 8.3 mm
 – Length 150 cm