



Laerdal Silicone Resuscitator (LSR)

The Laerdal Silicone Resuscitator (LSR) enables medical personnel to facilitate assisted ventilation of patients in pre-hospital, in-hospital and post-hospital care. The LSR comes in adult, paediatric and preterm sizes.

Adult

- 870050xx LSR Adult Basic without Mask
in Carton
- 870051xx LSR Adult Complete with Masks,
in Carton
- 870052xx LSR Adult Standard with Mask,
in Carton
- 870053xx LSR Adult Complete with Masks,
in Compact Case
- 870055xx LSR Adult Complete with Masks,
in Display Case

Paediatric

- 860050xx LSR Paediatric Basic without Mask,
in Carton
- 860051xx LSR Paediatric Complete with Masks,
in Carton

- 860052xx LSR Paediatric Standard
Child with Mask, in Carton
- 860053xx LSR Paediatric Complete with Masks,
in Compact Case
- 860055xx LSR Paediatric Complete with Masks,
in Display Case
- 860056xx LSR Paediatric Standard.Term with Mask,
in Carton

Preterm

- 850050xx LSR Preterm Basic without Mask,
in Carton
- 850051xx LSR Preterm Complete with Masks,
in Carton
- 850053xx LSR Preterm Complete with Masks
in Compact Case
- 850055xx LSR Preterm Complete with Masks
in Display Case

xx = language code

Models, Contents, and Containers

	ADULT (Patients >55 lbs.)			PEDIATRIC (Patients 5.5 to 55lbs.)				PRETERM (Patients < 5.5 lbs.)	
	Basic w/o Mask	Standard	Complete	Basic w/o Mask	Standard Child	Standard Term-Baby	Complete	Basic w/o Mask	Complete
Carton	870050	870052	870051	860050	860052	860056	860051	850050	850051
Compact Case			870053				860053		850053
Display Case			870055				860055		850055
Bag 1600 ml	*	*	*						
500 ml				*	*	*	*		
240 ml								*	*
Pt Valve w/o Relief Valve	*	*	*						
w/ Relief Valve				*	*	*	*	*	*
Mask 0/0						*	*		*
0/1						*	*		*
2							*		
3/4 w/Mask Cover			*				*		
4/5+ w/Mask Cover		*	*						
Intake/Reservoir Valve	*	*	*	*	*	*	*	*	*
Reservoir Bag 2600	*	*	*						
600				*	*	*	*	*	*
DFU	*	*	*	*	*	*	*	*	*

Laerdal Silicone Masks

Depending on the shape and size of your patient, there is a Laerdal Silicone mask that fits. An airtight seal between the skin and the silicone enables the proper administration of air. The Silicone masks are completely transparent for visibility of the patients' nose and mouth, fit standard 15/22mm connectors and are compatible with the Multi-Function Mask Cover.



- 860220 Child Silicone Mask 3-4
with Multi-function Mask Cover
- 860221 Child Silicone Mask 3-4 without Mask Cover
- 865200 Multi Function Mask Cover for Mask 3-4
- 870220 Adult Silicone Mask 4-5+
with Multi-function Mask Cover
- 870221 Adult Silicone Mask 4-5+ without Mask Cover
- 872220 Adult 4-5+ & Child 3-4 Silicone Mask
with Multi-function Mask Covers
- 875200 Multi-Function Mask Cover for Mask 4-5+

