



**STERIKING®**

**Sterile Barrier Systems & Products  
for Hospitals, Dental Practices & Laundries**

**WIPAK**


*We care that you pack safely!*

## See-Through Peel Pouches & Rolls

### S - Flat Pouches

For packing individual instruments, smaller sets and other items.\*

| Art.code | 2P | 3P | Size (mm) | Pieces/case |
|----------|----|----|-----------|-------------|
| S 17     | •  |    | 50 X 200  | 3.600       |
| S 1      | •  | •  | 50 X 250  | 2.400       |
| S 24     | •  | •  | 75 X 150  | 2.400       |
| S 2      | •  | •  | 75 X 200  | 2.400       |
| S 23     | •  |    | 75 X 230  | 1.800       |
| S 4      | •  | •  | 75 X 270  | 3.600       |
| S 18     | •  | •  | 75 X 300  | 4.200       |
| S 22     | •  |    | 75 X 520  | 2.400       |
| S 25     | •  |    | 100 X 150 | 2.400       |
| S 3      | •  | •  | 100 X 200 | 1.800       |
| S 5      | •  | •  | 100 X 270 | 3.000       |
| S 8      | •  | •  | 100 X 300 | 1.200       |
| S 19     | •  |    | 100 X 350 | 1.200       |
| S 9      | •  | •  | 100 X 400 | 1.200       |
| S 12     | •  |    | 100 X 570 | 1.800       |
| S 27     | •  |    | 120 X 400 | 1.200       |
| S 15     | •  | •  | 150 X 200 | 1.200       |
| S 6      | •  | •  | 150 X 270 | 1.800       |
| S 20     | •  | •  | 150 X 300 | 1.800       |
| S 26     | •  | •  | 150 X 350 | 1.800       |
| S 13     | •  | •  | 150 X 400 | 1.200       |
| S 28     | •  |    | 150 X 520 | 1.200       |
| S 34     | •  |    | 160 X 600 | 1.200       |
| S 33     | •  |    | 180 X 320 | 1.000       |
| S 7      | •  | •  | 205 X 270 | 600         |
| S 7 A    | •  |    | 205 X 340 | 1000        |
| S 10     | •  | •  | 205 X 400 | 600         |
| S 10 A   | •  |    | 205 X 440 | 600         |
| S 35     | •  |    | 205 X 600 | 1.200       |
| S 2070   | •  |    | 205 X 700 | 600         |
| S 16     | •  | •  | 250 X 380 | 1.200       |
| S 11     | •  | •  | 250 X 500 | 600         |
| S 36     | •  |    | 250 X 600 | 500         |
| S 37     | •  |    | 250 X 700 | 600         |
| S 29     | •  |    | 270 X 350 | 600         |
| S 30     | •  | •  | 270 X 440 | 600         |
| S 21     | •  | •  | 300 X 500 | 600         |
| S 14     | •  | •  | 300 X 570 | 600         |
| S 38     | •  |    | 320 X 500 | 600         |
| S 31     | •  |    | 420 X 500 | 500         |
| S 32     | •  | •  | 420 X 600 | 500         |



### B - Pouches with Gusset

For more voluminous items.\*

| Art.code | 2P | Size (mm)       | Pieces/case |
|----------|----|-----------------|-------------|
| B 35     | •  | 100 X 50 X 300  | 500         |
| B 30     | •  | 100 X 50 X 360  | 500         |
| B 31     | •  | 150 X 50 X 400  | 500         |
| B 32     | •  | 150 X 50 X 460  | 500         |
| B 36     | •  | 200 X 55 X 400  | 250         |
| B 33     | •  | 200 X 55 X 500  | 250         |
| B 37     | •  | 250 X 65 X 480  | 250         |
| B 34     | •  | 300 X 80 X 550  | 250         |
| B 38     | •  | 420 X 100 X 570 | 250         |

### R - Flat Rolls \*

| Art.code | 2P | 3P | Size (mm)   | Pieces/case |
|----------|----|----|-------------|-------------|
| R 39     | •  | •  | 50 : 200 m  | 5           |
| R 40     | •  | •  | 75 : 200 m  | 4           |
| R 41     | •  | •  | 100 : 200 m | 3           |
| R 41A    | •  | •  | 120 : 200 m | 2           |
| R 42     | •  | •  | 150 : 200 m | 2           |
| R 43     | •  | •  | 200 : 200 m | 1           |
| R 44     | •  | •  | 250 : 200 m | 1           |
| R 45     | •  | •  | 300 : 200 m | 1           |
| R 46     | •  | •  | 350 : 200 m | 1           |
| R 47     | •  | •  | 400 : 200 m | 1           |

| Art.code | 2P | 3P | Size (mm)   | Pieces/case |
|----------|----|----|-------------|-------------|
| R 39/100 | •  |    | 50 : 100 m  | 8           |
| R 40/100 | •  |    | 75 : 100 m  | 8           |
| R 41/100 | •  |    | 100 : 100 m | 4           |
| R 42/100 | •  |    | 150 : 100 m | 2           |
| R 43/100 | •  |    | 200 : 100 m | 2           |
| R 44/100 | •  |    | 250 : 100 m | 2           |
| R 45/100 | •  |    | 300 : 100 m | 1           |
| R 47/100 | •  |    | 400 : 100 m | 1           |
| R 49/100 | •  |    | 500 : 100 m | 2           |

| Art.code | 2P | 3P | Size (mm)  | Pieces/case |
|----------|----|----|------------|-------------|
| R 39/31  | •  |    | 50 : 31 m  | 10          |
| R 40/31  | •  |    | 75 : 31 m  | 10          |
| R 41/31  | •  |    | 100 : 31 m | 10          |
| R 42/31  | •  |    | 150 : 31 m | 10          |
| R 43/31  | •  |    | 200 : 31 m | 6           |
| R 44/31  | •  |    | 250 : 31 m | 6           |
| R 45/31  | •  |    | 300 : 31 m | 6           |

### RB - Rolls with Gusset \*

| Art.code | 2P | 3P | Size (mm)        | Pieces/case |
|----------|----|----|------------------|-------------|
| RB 50    | •  | •  | 75 X 25 : 100 m  | 4           |
| RB 51    | •  | •  | 100 X 50 : 100 m | 3           |
| RB 52    | •  | •  | 150 X 50 : 100 m | 2           |
| RB 53    | •  | •  | 200 X 55 : 100 m | 1           |
| RB 54    | •  | •  | 250 X 65 : 100 m | 1           |
| RB 55    | •  | •  | 300 X 80 : 100 m | 1           |
| RB 56    | •  | •  | 350 X 80 : 100 m | 1           |
| RB 57    | •  | •  | 400 X 80 : 100 m | 1           |

\* 2P = steam and EO gas, 3P = steam, EO and FO.

## SMX - High Performance ProWraps

### SMX-1 Blue Light

For small and medium-size trays and sets.

| Art.code | Size (mm)   | Pieces/case |
|----------|-------------|-------------|
| SMX 1060 | 600 x 600   | 480         |
| SMX 1075 | 750 x 750   | 200         |
| SMX 1090 | 900 x 900   | 160         |
| SMX 1100 | 1000 x 1000 | 160         |
| SMX 1120 | 1200 x 1200 | 120         |

### SMX-2 Blue Regular

For medium and large-size trays and sets and for container inner wraps.

| Art.code | Size (mm)   | Pieces/case |
|----------|-------------|-------------|
| SMX 2060 | 600 x 600   | 480         |
| SMX 2075 | 750 x 750   | 200         |
| SMX 2090 | 900 x 900   | 120         |
| SMX 2100 | 1000 x 1000 | 120         |
| SMX 2120 | 1200 x 1200 | 120         |

### SMX-2+ Green Regular

Optimal choice for inner wrap, as it offers a reliable sterile field.

| Art.code  | Size (mm)   | Pieces/case |
|-----------|-------------|-------------|
| SMX 2+060 | 600 x 600   | 480         |
| SMX 2+075 | 750 x 750   | 200         |
| SMX 2+090 | 900 x 900   | 120         |
| SMX 2+100 | 1000 x 1000 | 120         |
| SMX 2+120 | 1200 x 1200 | 120         |

### SMX-3 Blue Heavy Duty

Heavier structure for the most demanding applications.

| Art.code | Size (mm)   | Pieces/case |
|----------|-------------|-------------|
| SMX 3075 | 750 x 750   | 120         |
| SMX 3090 | 900 x 900   | 120         |
| SMX 3100 | 1000 x 1000 | 120         |
| SMX 3120 | 1200 x 1200 | 104         |

### SMX-4 Blue Super Duty

Especially high thickness for orthopedic sets and other heavy-duty uses.

| Art.code | Size (mm)   | Pieces/case |
|----------|-------------|-------------|
| SMX 4075 | 750 x 750   | 120         |
| SMX 4090 | 900 x 900   | 104         |
| SMX 4100 | 1000 x 1000 | 104         |
| SMX 4120 | 1200 x 1200 | 80          |

## NWB - Nonwoven Soft Blue

For general use, for laundry and instrument sets.

| Art.code | Size (mm)   | Pieces/case |
|----------|-------------|-------------|
| NWB 50   | 500 x 500   | 500         |
| NWB 60   | 600 x 600   | 500         |
| NWB 75   | 750 x 750   | 250         |
| NWB 90   | 900 x 900   | 250         |
| NWB 100  | 1000 x 1000 | 250         |
| NWB 120  | 1200 x 1200 | 100         |
| NWB 137  | 1370 x 1370 | 100         |

## NWG - Nonwoven Super Soft Green

Superior softness and toughness for demanding applications.

| Art.code                                     | Size (mm)            | Pieces/case |
|--|----------------------|-------------|
| NWG 50                                       | 500 x 500            | 400         |
| NWG 60                                       | 600 x 600            | 400         |
| NWG 75                                       | 750 x 750            | 200         |
| NWG 90                                       | 900 x 900            | 150         |
| NWG 100                                      | 1000 x 1000          | 150         |
| NWG 120                                      | 1200 x 1200          | 100         |
| Available as interleaved with SMX Wrap-grade |                      |             |
| NWGL   | Nonwoven Green Light |             |



## SPC - Sterilization Paper Crepe

Special grade paper creped to give increased flexibility.

| Art.code    | Size (mm)   | Pieces/case |
|-------------|-------------|-------------|
| SPC 40      | 400 x 400   | 500         |
| SPC 50      | 500 x 500   | 500         |
| SPC 60      | 600 x 600   | 500         |
| SPC 75      | 750 x 750   | 250         |
| SPC 90      | 900 x 900   | 250         |
| SPC 100     | 1000 x 1000 | 250         |
| SPC 120     | 1200 x 1200 | 100         |
| SPC 100/140 | 1000 x 1400 | 100         |

## SCB - Super Crepe Blue

Soft crepe - for improved drapeability.

| Art.code | Size (mm)   | Pieces/case |
|----------|-------------|-------------|
| SCB 50   | 500 x 500   | 500         |
| SCB 60   | 600 x 600   | 500         |
| SCB 75   | 750 x 750   | 250         |
| SCB 90   | 900 x 900   | 250         |
| SCB 100  | 1000 x 1000 | 250         |
| SCB 120  | 1200 x 1200 | 100         |

## Absorbent Trayliners

For use in surgical instrument trays to absorb excessive condensation.

| TL-Standard               | Size (mm) | Pieces/case |
|---------------------------|-----------|-------------|
| TL 25                     | 250 x 300 | 1000        |
| TL 40                     | 400 x 500 | 500         |
| TL 50                     | 500 x 300 | 500         |
| TLWE-Trayliners with edge |           |             |
| TLWE 2525                 | 250 x 250 | 500         |
| TLWE 4825                 | 480 x 250 | 500         |
| TLWE 5030                 | 500 x 300 | 500         |
| TLWE-Trayliners nonwoven  |           |             |
| TLNW 1525                 | 150 x 250 | 3.200       |
| TLNW 3025                 | 300 x 250 | 1.200       |
| TLNW 3040                 | 300 x 400 | 1.200       |
| TLNW 3050                 | 300 x 500 | 600         |



## Indicator Strips

AC = steam, EO-gas, FO-gas and  
PL = low temperature plasma.

NEW!

| Art.code   | Strips/dispenser | Pieces/case |
|------------|------------------|-------------|
| AC 6       | 250              |             |
| AC 250/500 | 250              | 1.000       |
| EO 250     | 250              | 1.000       |
| FO 250/500 | 250              |             |
| PL 500     |                  | 500         |

## Daily Test Sheet

For Bowie & Dick Test with a cloth package.

| Art.code | Pieces/case |
|----------|-------------|
| DT 50    | 50          |

## Daily Test Pack

The ready made test packs for Bowie & Dick Test.

| Art.code | Pieces/case |
|----------|-------------|
| DCC2     | 20          |
| DC 20    | 20          |

## Helix Challenge Test

For a charge control and as a steam penetration test for hollow instruments with small lumina.

| Art.code | Strips/dispenser | Pieces/case |
|----------|------------------|-------------|
| HX 100   | 100              | 25          |
| HX 400   | 400              | 25          |

## Tapes

Indicator Tapes, S=steam, E=EO-gas and H=hot-air indicator

| Art.code              | Width/length | Pieces/case |
|-----------------------|--------------|-------------|
| ITS 12 G              | 12 mm x 50 m | 72          |
| ITS 19 G              | 19 mm x 50 m | 48          |
| ITS 25 G              | 25 mm x 50 m | 36          |
| JTB 19 S G (B = blue) | 19 mm x 50 m | 48          |
| JTB 25 S G (B = blue) | 25 mm x 50 m | 36          |
| ITE 19                | 19 mm x 50 m | 48          |
| ITH 19                | 19 mm x 50 m | 48          |

Closing Tapes with no imprint

| Art.code            | Width/length | Pieces/case |
|---------------------|--------------|-------------|
| ITN 19              | 19 mm x 50 m | 48          |
| ITN 25              | 25 mm x 50 m | 36          |
| JTB 19 N (B = blue) | 19 mm x 50 m | 48          |
| JTB 25 N (B = blue) | 25 mm x 50 m | 36          |

Plasma Indicator Tape Strip of DuPont™ Tyvek®

| Art.code | Width/length   | Pieces/case |
|----------|----------------|-------------|
| ITP 20   | 20 mm x 150 mm | 200         |

Tape dispenser for one roll

| Art.code | Pieces/case |
|----------|-------------|
| TDS 25   | 1           |

## SS - Standard Self Seal Pouches

with adhesive closure. 2P = steam and EO gas. **NEW UNITS!**

| Art.code | Size (mm) | Pieces/case | Pieces/case |
|----------|-----------|-------------|-------------|
| SS 1     | 90 x 200  | 1.200       | 1.000       |
| SS 2     | 90 x 250  | 1.200       | 1.000       |
| SS 3     | 90 x 570  | 1.200       | 1.000       |
| SS 4     | 130 x 270 | 1.200       | 1.000       |
| SS 4 A   | 130 x 380 | 1.200       | 1.000       |
| SS 5     | 200 x 350 | 1.200       | 1.200       |
| SS 5 A   | 190 x 330 | 1.200       | 1000        |
| SS 6     | 250 x 400 | 1.200       | 600         |
| SS 7     | 300 x 450 | 600         | 600         |

## SST - Self Seal Plus Pouches

with an attached internal indicator strip.

| Art.code | Size (mm) | Pieces/case |
|----------|-----------|-------------|
| SST 1    | 90 x 200  | 1.000       |
| SST 2    | 90 x 200  | 1.000       |
| SST 3A   | 100 x 640 | 1.000       |
| SST 4    | 130 x 270 | 1.000       |
| SST 4A   | 130 x 380 | 1.000       |
| SST 5A   | 190 x 330 | 1.000       |
| SST 6    | 250 x 400 | 600         |
| SST 7    | 300 x 450 | 600         |

## SNW - TEX-Line Pouches

Peel pouches with nonwoven and film for packing heavy and bulky items. 3P = steam, EO and FO

| Art.code | Size (mm) | Pieces/case |
|----------|-----------|-------------|
| SNW 2040 | 205 x 400 | 500         |
| SNW 2744 | 270 x 440 | 500         |
| SNW 3250 | 320 x 500 | 500         |
| SNW 4260 | 420 x 600 | 300         |
| SNW 5065 | 500 x 650 | 200         |
| SNW 5075 | 500 x 750 | 200         |
| SNW 5775 | 570 x 750 | 200         |

## CB - Dust Cover Bags

for use after sterilization.

| Art.code | Size (mm) | Pieces/case |
|----------|-----------|-------------|
| CB 1     | 150 x 270 | 1.000       |
| CB 2     | 205 x 300 | 1.000       |
| CB 3     | 250 x 400 | 1.000       |
| CB 4     | 300 x 500 | 500         |
| CB 5     | 400 x 550 | 500         |
| CB 6     | 420 x 700 | 500         |

## SSDC - Self Seal Dust Cover Bags

with Adhesive Closure

| Art.code | Size (mm) | Pieces/case |
|----------|-----------|-------------|
| SSDC 10  | 150 x 270 | 500         |
| SSDC 20  | 200 x 350 | 500         |
| SSDC 30  | 250 x 400 | 500         |
| SSDC 40  | 300 x 550 | 200         |
| SSDC 60  | 420 x 700 | 200         |

## LTS - Flat Tyvek® Pouches

Made with DuPont™ Tyvek® for low temperature sterilization.

Printed with hydrogen plasma process indicator. **NEW UNITS!**


| Art.code | Size (mm) | Pieces/case | Pieces/case |
|----------|-----------|-------------|-------------|
| LTS 7520 | 75 x 200  | 1.200       | 1.000       |
| LTS 1025 | 100 x 250 | 1.200       | 1.000       |
| LTS 1530 | 150 x 300 | 1.200       | 1.000       |
| LTS 1644 | 160 x 440 | 1.200       | 600         |
| LTS 1660 | 160 x 600 | 600         | 600         |
| LTS 2038 | 205 x 380 | 1.200       | 1.200       |
| LTS 2538 | 250 x 380 | 1.200       | 600         |
| LTS 2550 | 250 x 500 | 600         | 600         |

## LTSS - Self Seal Tyvek® Pouches

| Art.code | Size (mm) | Pieces/case |
|----------|-----------|-------------|
| LTSS 1   | 90 x 200  | 1.000       |
| LTSS 2   | 90 x 250  | 1.000       |
| LTSS 4   | 130 x 270 | 1.000       |
| LTSS 4 A | 130 x 380 | 1.000       |
| LTSS 5 A | 190 x 330 | 1.000       |
| LTSS 6   | 250 x 400 | 600         |
| LTSS 7   | 300 x 450 | 600         |

## LTR - Flat Tyvek® Rolls

| Art.code | Size (mm) | Pieces/case |
|----------|-----------|-------------|
| LTR 40   | 75 x 100  | 1           |
| LTR 41   | 100 x 100 | 1           |
| LTR 42   | 150 x 100 | 1           |
| LTR 43   | 200 x 100 | 1           |
| LTR 43 A | 225 x 100 | 1           |
| LTR 44   | 250 x 100 | 1           |
| LTR 45   | 300 x 100 | 1           |
| LTR 46   | 350 x 100 | 1           |
| LTR 47   | 400 x 100 | 1           |
| LTR 49   | 500 x 100 | 1           |





## Sealing Machines

**NEW!** Model: RS 120A Rotosealer

Model: RS 220A Advanced Rotosealer with printer

Model: IS 250 Impulse sealer with cutting device

Cutting Devices in 3 different widths

• CD 600, CD 900, CD 1200

## Seal Control Sheet

For operational qualification of the sealing process.

| Art.code        | Strips/dispenser | Pieces/case |
|-----------------|------------------|-------------|
| SC 250 (paper)  | 250              | 1.000       |
| TSC 200 (Tyvek) | 200              | 800         |

**NEW!**

## Non-toxic Markers

Waterproof, heat-resistant markers for sterilization packages. 10/dispenser.

| Art.code | Color Pieces | Size (mm) | Pieces/case |
|----------|--------------|-----------|-------------|
| MAR 1    | ■            | 0.75      | 100         |
| MAR 1T   | ■            | 1.00      | 100         |
| MAR 2    | ■            | 0.75      | 100         |
| MAR 3    | ■            | 0.75      | 100         |
| MAR 4    | ■            | 0.75      | 100         |

## PB - Paper Bags with Gusset

Heat-sealable lacquer top and steam indicator.

| Art.code | Size (mm)       | Pieces/case |
|----------|-----------------|-------------|
| PB 1     | 90 x 50 x 125   | 1.000       |
| PB 2     | 110 x 30 x 190  | 1.000       |
| PB 3     | 125 x 50 x 250  | 1.000       |
| PB 4     | 140 x 75 x 250  | 500         |
| PB 5     | 140 x 50 x 330  | 500         |
| PB 6     | 165 x 50 x 380  | 500         |
| PB 7     | 180 x 95 x 380  | 500         |
| PB 8     | 190 x 65 x 330  | 500         |
| PB 9     | 250 x 100 x 380 | 500         |
| PB 10    | 300 x 75 x 530  | 250         |

## HR - Hot Air Rolls

Special polyamide tubing for sterilization by dry heat.

| Art.code | Width/length   | Pieces/case |
|----------|----------------|-------------|
| HR 39    | 50 mm x 200 m  | 5           |
| HR 40    | 75 mm x 200 m  | 4           |
| HR 41    | 100 mm x 200 m | 3           |
| HR 42    | 150 mm x 200 m | 2           |
| HR 43    | 200 mm x 200 m | 1           |
| HR 44    | 250 mm x 200 m | 1           |

## Registration Labels

For registering the sterilized items, enabling recall of all items if necessary.

PP-based labels, resistant to disinfection, washing and sterilization

| Art.code |                      | Pieces/case |
|----------|----------------------|-------------|
| TRLAB    | Tray Label unprinted | 1.000       |
| TRLABPR  | Tray Label printed   | 1.000       |

Paper label with adhesive backing, suitable as a bar code carrier

| Art.code |           | Pieces/case |
|----------|-----------|-------------|
| TRTAG    | Tray Tag  | 1.000       |
| ENV 1    | Envelopes | 1.000       |

# STERIKING®

Wipak Oy  
P.O. Box 45  
Wipaktie 2  
FI-15561 Nastola, Finland  
Tel: +358-20510 311  
Fax: +358-20510 3333  
www.wipak.com  
www.steriking.info  
e-Mail: steriking@wipak.com



# WIPAK

We reserve the right to make changes in the products.