

# Steriliseringspåsar - med bälg

## Steriliseringspåse papier/plast bälg



Art.nr.	Storlek	Förp.
WI-B35	100 x 50 x 300 mm	50 / 500
WI-B30	100 x 50 x 360 mm	50 / 500
WI-B31	150 x 50 x 400 mm	50 / 500
WI-B32	150 x 50 x 460 mm	50 / 500
WI-B36	200 x 55 x 400 mm	50 / 250
WI-B33	200 x 55 x 500 mm	50 / 250
WI-B37	250 x 65 x 480 mm	50 / 250
WI-B34	300 x 80 x 550 mm	50 / 250
WI-B38	420 x 100 x 570 mm	50 / 250

## Steriliseringspåse papier/plast, bälg på rulle



Art.nr.	Storlek	Förp.
WI-RB50	75 x 20 mm x 100 m	1 / 4
WI-RB51	100 x 50 mm x 100 m	1 / 3
WI-RB52	150 x 50 mm x 100 m	1 / 2
WI-RB53	200 x 55 mm x 100 m	1 / 1
WI-RB54	250 x 65 mm x 100 m	1 / 1
WI-RB55	300 x 80 mm x 100 m	1 / 1
WI-RB56	350 x 80 mm x 100 m	1 / 1
WI-RB57	400 x 80 mm x 100 m	1 / 1