



STERIKING® Sterimarkers

STERIKING® Sterimarkers are developed for permanent reliable marking of sterilization packages on the film side, on tape or on a tray label. STERIKING® Sterimarkers have a special ink (available in black, red, blue and green) which is immediately dry, waterproof, light proof, resistant to atmospheric action, heat resistant and acid resistant. The ink is not harmful for the intended use and can be sterilized in steam, EO, FO, VH₂O₂ and by irradiation.

Storage Recommendations

It is recommended that the STERIKING® products are kept in the original, closed transport carton and are stored in dry and clean conditions protected from direct sunlight and excessive moisture.

Do not shake the markers and keep the cap on the marker when not in use.

Sales and Transport Packing

STERIKING® Sterimarkers are packed in white boxes, 10 pieces per box.

Please refer to national and local regulations regarding waste disposal.

In Case of Complaint

In event of any complaint, the lot number and identification code must be provided by the complainant. For evaluation of claimed product, a defective sample (or a digital photo) and description of the defect together with an unused specimen must be made available to Wipak.

STERIKING® is a registered trademark of Wipak Medical.

Steriking® Sterimarkers

Code	Colour	Thickness	Markers/box
MAR1	Black	0,75	10
MAR2	Red	0,75	10
MAR3	Blue	0,75	10
MAR4	Green	0,75	10
MAR1T	Black	1,0	10

Min. order quantity 10 boxes

This specification refers to the named product group and shall be valid until the next revision. Other product related documents may be available upon request.

The information contained here is to our knowledge accurate and reliable as of the date of the publication. Wipak extends no warranties and makes no representations as to the accuracy or completeness of the information contained herein, and assumes no responsibility regarding the consequences of its use or for any printing errors. It is the customer's responsibility to inspect and test our products in order to satisfy himself as to the suitability of the products for the customer's particular purpose and suitability to the actual circumstances the product is exposed to. The customer is also responsible for the appropriate, safe and legal use, processing and handling of our products especially when recommendations for safe use and storage are given. Nothing herein shall constitute any warranty, nor is a protection from any law or patent to be inferred.