

BD MaxPlus™ Clear needle-free connector

SKU/REF MP1000C-0006
SKU/REF Name BD MaxPlus™ Clear needle-free connector



Specifications

Length (cm)	3.3
Priming Volume (mL)	0.28
Maximum Pressure Rating	325 PSI@10 mL per second
Needle-Free Connector	BD MaxPlus™
Change Interval	7 days/200 activations
DEHP	DEHP-Free
Sterilization Method	Electron Beam
Country of Manufacture	Mexico
Classification	Ila
Quantity per Box	100
Quantity per Case	2400
Shelf Life	3 Years

CE 2797

<https://www.bd.com/en-eu/bd-maxplus-clear-needle-free-connector/mp1000c-0006?return=true>

Becton Dickinson Sverige AB, Marieviksgatan 25, 11743 Stockholm
bdsweden@bd.com

bd.com

BD, BD-logo, NeutroX, NeutraClear är varumärken som tillhör Becton, Dickinson and Company eller dess dotterbolag
© 2022 BD. Alla rättigheter förbehållna. BD-59783 March 2022

