

Giving breath to life

The Vital Signs® T-Piece

Resuscitation can be crucial to a newborn's wellness in the moments after birth. The Vital Signs T-Piece resuscitator kit is a single-patient-use resuscitator circuit intended to provide respiratory support in such situations. They can be used in conjunction with equipment, such as the GE Giraffe™ warmers, or a stand-alone infant resuscitation system.

All T-piece circuits include a positive end-expiratory pressure (PEEP) valve to help maintain optimal pressure. Made of corrugated tubing without natural rubber latex or DEHP, the T-Piece circuit kits can include a size 0 face mask, size 1 face mask or no mask. Masks feature air-filled cushions and an anatomical shape with finger ports. A dust cap with a strap is also included.



Product components

M1091335VS	T-piece neonatal patient circuit kit, disposable, 10 per package
M1091316VS	T-piece neonatal patient circuit kit, disposable, with size 0 mask, 10 per package
M1091365VS	T-piece neonatal patient circuit kit, disposable, with size 1 mask, 10 per package

Vital Signs®

TIMIK
MEDICAL

Product specifications

Tube length	183 cm (193 cm with T-piece connector)
Tube diameter	10 mm
Dead space	<5 ml
Patient connection	ISO 5356-1 15 mm taper, female

Single-patient use

Material	Product does not contain natural rubber latex or DEHP
----------	---

Intended for use with neonates and infants weighing less than 10 kg (22 lb)

Individually packed circuits

Product weight

M1091335VS	46 g (ea.)
M1091316VS	58 g (ea.)
M1091365VS	68 g (ea.)

