

StatLock® Universal Plus Stabilisation Devices

- Minimise catheter motion and reduce overall complications
- Improve patient comfort by eliminating suture securement
- Prevent needlestick injuries associated with catheter securement³

Provides solutions that follow OSHA's recommendations for sutureless securement³



VUPD68

Device Retainers

- Releasable device clamp
- When closed, allows for securement of lubricious catheters
- When open, permits gradual withdrawal of catheter without pad replacement
- Four sizes offers versatility by fitting most catheter sizes

Device Anchor Pads

- Custom device pad (with underlying hydrocolloid swath) allows for skin-safe adherence
- Allows for visibility of insertion site
- Alcohol-soluble adhesive is gentle on the skin
- Breathable, latex-free Tricot fabric provides added catheter security and patient comfort
- Includes skin protectant for enhanced stabilisation

Description	Product Code	Box Quantity
StatLock Universal Plus, Butterfly Tricot Pad, Small (6-8.5Fr.)	VUPD68	25
StatLock Universal Plus, Butterfly Tricot Pad, Medium (10-12Fr.)	VUPD1012	25
StatLock Universal Plus, Butterfly Tricot Pad, Large (12-14Fr.)	VUPD1214	25
StatLock Universal Plus, Butterfly Tricot Pad, Extra Large (14-16Fr.)	VUPD1416	25