

# AirLife® neonatal care

Helping meet the challenging clinical needs of our most critical newborns

Each newborn has a unique set of needs that must be assessed quickly and efficiently. Our innovative AirLife® neonatal intensive care unit (NICU) supplies and accessories—designed for use with GE Healthcare incubators, warmers, resuscitation systems, and BiliSoft™ phototherapy systems—help enable clinicians to provide exceptional care to the most delicate neonatal patients.



**AirLife®**



### T-piece resuscitator kits

Resuscitation can be crucial to a newborn's health in the moments after birth. The T-piece resuscitator kit is a single-patient-use resuscitator circuit that provides respiratory support.

All T-piece circuits include a positive end-expiratory pressure (PEEP) valve to help maintain optimal pressure. Made of corrugated tubing without latex or added DEHP, the T-piece circuit kits come with a choice of face mask (*either size 0 or 1*) or no mask. The attached dust cap helps to minimize risk of contamination.

Product no.	Description	Qty. (cs.)
M1091335VS	T-piece neonatal circuit	10
M1091316VS	T-piece neonatal circuit kit, size 0 mask	10
M1091365VS	T-piece neonatal circuit kit, size 1 mask	10



### Infant skin temperature probes

Temperature is a key vital sign when monitoring premature infants, as they are vulnerable to having abnormally low core temperatures. Our infant skin temperature probes provide a less invasive method of temperature management to help monitor the body temperature of premature newborns who experience difficulty maintaining normothermia.

Product no.	Description	Qty. (cs.)
2074816-001	Infant skin temperature probe, disposable	10
2074817-001	Infant skin temperature probe, disposable	50
2075796-001	Infant skin temperature probe, reusable	1 ea.
0203-1980-300	Heat reflective patch, for use with reusable infant temp probe	50



## Phototherapy covers

AirLife disposable phototherapy covers help provide simple intensive therapy when newborns present with hyperbilirubinemia. These ultra-soft disposable pad covers and nests are made of flame-retardant fabric gentle on sensitive newborn skin. Designed to help promote skin to skin contact, the pad covers and nests enable comfortable swaddling and much needed support during phototherapy treatments.

Product no.	Description	Qty. (cs.)
For babies less than 1,500 g		
M1093120	Small disposable pad cover	50
M1097108	Small disposable pad cover	20
M1093122	Small disposable nest	15
For babies at least 1,500 g or more		
M1093121	Large disposable pad cover	50
M1097109	Large disposable pad cover	20
M1093123	Large disposable nest	15

Disposable nest offers boundaries and support with a cushy foot roll and gentle, transparent straps.



Disposable pad cover slips over the LED pad with transparent straps to wrap around the baby while swaddling.





To learn more, contact your CareFusion sales representative,  
call **1.800.932.0760, option 2** or visit [carefusion.com/AirLife](http://carefusion.com/AirLife).

CareFusion  
Vernon Hills, IL

[carefusion.com](http://carefusion.com)

