




**FRESENIUS  
KABI**

**Technical Customer Specification**

Compoflex® 1F 63 ml CPDA1 - PDS-V

Nr.:	PS11150/ 7
Language(s):	none
Page	3 of 3
Released:	17-5-2016
Shelf life:	24 months



<b>Label:</b>	PVC / Barcode according ISBT128B	<b>Plastic Foil:</b>	PVC - DEHP	<b>Container Volume:</b>	600 ml	<b>Donation Volume:</b>	450 ml		
<b>Manufacturer:</b>	Fresenius Kabi AG, 61346 Bad Homburg, Germany	<b>Needle:</b>	16G - Trifacial Bevel / Retractable needle with optimized flow						
<b>Sterilization:</b>	Steam sterilization	<b>Primary:</b>	1 unit	<b>Secondary:</b>	9 units	<b>Carton:</b>	36 units	<b>Max. Pallet:</b>	1080 units