






Double swivel elbows





Code	1897	1899	1898
Box qty	35	40	35
Flip top cap		✓	
7.6mm port		✓	
Double flip top cap with seal			✓
Sterile option code	1897000S	1899000S	1898000S
Sterile box qty	15	15	15

Fixed elbows for direct patient connection



Code	1997	1993	2714	1883	1992	1995
Box qty	35	35	25	70	25	25
Luer lock port			✓	✓		
Flip top cap		✓				
7.6mm port		✓				✓
Sterile option code	1997000S		2714000S			
Sterile box qty	15		15			



Code	1892	1893	1891	1990	1991
Box qty	25	25	25	25	20
7.6mm port	✓	✓	✓		
6mm stem				✓	✓