










PRODUCT NAME	MAIN ARTICLE NUMBER	USE	CLASS/ STANDARD	UNIT OF MEASURE
WASH MONITORING				
Label Color		Test Soil		
 Test Soil	504229300	W-D	EN ISO 15883	10
Wash Monitoring and Holders				
 Wash Monitor	503864500	W-D	EN ISO 15883	25
 Wash Monitor	503864600	W-D	EN ISO 15883	50
 Wash Monitor	503864700	W-D	EN ISO 15883	100
 Wash Monitor FR • Flexible & Rigid	503864900	W-D	EN ISO 15883	25
 Wash Monitor Holder	503864800	W-D/Ultrasonic	EN ISO 15883	6
 Wash Monitor Holder R • Rigid Lumen Test	503865000	W-D/Ultrasonic	EN ISO 15883	1
 Wash Monitor Holder F • Flexible Lumen Test	503865100	W-D/Ultrasonic	EN ISO 15883	1
Wash Monitoring, Ultrasonic				
 Wash Monitor U • Ultrasonic Test	503865200	Ultrasonic	EN ISO 15883	50
 Wash Monitor U/L • Ultrasonic/Lumen Test	503865300	Ultrasonic	EN ISO 15883	25
Protein Tests				
 Protein Test Instrument Surface	503878500	W-D/Ultrasonic	EN ISO 15883	25
 Protein Test Instrument Lumen	503911200	W-D/Ultrasonic	EN ISO 15883	25
 Protein Test Flexible Endoscope 2,5m	503964400	W-D/Ultrasonic	EN ISO 15883	10
STEAM STERILIZATION PRODUCTS - DAILY MACHINE RELEASE				
Bowie-Dick Tests				
 Bowie-Dick Test Pack 134°C, Getinge Turbo Sterilizers	504052100	Steam	EN ISO 11140-4	30
 Getinge Check Bowie-Dick Test Pack 134°C, Sub-trans cycles	601154901	Steam	EN ISO 11140-4	30
 Bowie-Dick Test Pack 121°C	504052200	Steam	EN ISO 11140-4	30
 Bowie-Dick Test Card TT (Table Top Sterilizers)	504052300	Steam	BS EN 867-5	15
 Bowie-Dick Test Card Kit 134°C	504052500	Steam	EN ISO 11140-4	1 holder+31 tests
 Bowie-Dick Test Card 134°C	504052400	Steam	EN ISO 11140-4	120
 AdvaTec Bowie-Dick Test Pack, Getinge Turbo Sterilizers	504225100	Steam	EN ISO 11140-4	30
 Getinge Check Bowie-Dick Mini-Pack, Sub-trans cycles	6001155600	Steam	ISO 11138	30
STEAM STERILIZATION PRODUCTS - PACK RELEASE				
Emulating Indicators Non-Adhesive				
 Emulating Indicator 4	504051000	Steam	EN ISO 11140-1 Class 6	250
 Emulating Indicator 5.3	504050800	Steam	EN ISO 11140-1 Class 6	250
 Emulating Indicator 7	504050900	Steam	EN ISO 11140-1 Class 6	250
 Emulating Indicator 18	504051400	Steam	EN ISO 11140-1 Class 6	250
Emulating Indicators - Adhesive				
 Emulating Indicator 3.5 -A	504051100	Steam	EN ISO 11140-1 Class 6	250
 Emulating Indicator 4 -A	504052600	Steam	EN ISO 11140-1 Class 6	250
 Emulating Indicator 5.3 -A	504051200	Steam	EN ISO 11140-1 Class 6	250
 Emulating Indicator 7 -A	504051300	Steam	EN ISO 11140-1 Class 6	250
 Multi-Parameter Indicator	504051700	Steam	EN ISO 11140-1 Class 4	250
 Multi-Parameter Indicator L	504051800	Steam	EN ISO 11140-1 Class 4	250
STEAM STERILIZATION PRODUCTS - LOAD RELEASE				
 Helix Test 3.5	504052700	Steam	BS EN 867-5 Class 2	1 device + 250 tests
 Helix Test 4	504052800	Steam	Class 2	1 device + 250 tests
 Helix Test 5.3	504052900	Steam	Class 2	1 device + 250 tests
 Helix Test 7	504053000	Steam	Class 2	1 device + 250 tests
 Tape 18 mm	504053200	Steam	EN ISO 11140-1 Class 1	48
 Tape 24 mm	504053300	Steam	EN ISO 11140-1 Class 1	36
 Green Tape 18 mm	6001154902	Steam	EN ISO 11140-1 Class 1	48
 Green Tape 24 mm	6001154903	Steam	EN ISO 11140-1 Class 1	36
Biological Indicators				
 Self-Contained Biological Indicator	504053900	Steam	ISO 11138	100
 Self-Contained Biological Indicator	504054000	EO	ISO 11138	100
 Smart-Read Biological Test Pack	504054500	Steam	ANSI/AAMI ST79, ISO 11138, EN ISO 11140-1 Class 6	30 tests + 5 controls
 Multi-Incubator	504054700	24 hrs	-	1
 Smart-Read Biological Indicator	504054100	Steam	ISO 11138	100
 Self-Contained Biological Test Pack	504054600	Steam	ANSI/AAMI ST79, ISO 11138, EN ISO 11140-1 Class 6	30 tests + 5 controls
 Smart Well Incubator with Printer	504054400	10 hrs	-	1
SPECIALTY PRODUCTS				
 Helix Test	504053100	FORM	-	1 device + 250 tests
 Multi-Parameter Indicator	504051900	FORM	EN ISO 11140-1 Class 4	250
 Tape 18 mm	504053600	EO	EN ISO 11140-1 Class 1	48
 Tape 24 mm	504053700	EO	EN ISO 11140-1 Class 1	36
 Multi-Parameter Indicator	504051600	EO	EN ISO 11140-1 Class 4	250
 Multi-Parameter Indicator	504054800	Dry Heat	EN ISO 11140-1 Class 4	250
 Plasma Indicator	504052000	Plasma	-	250