

PRODUKTINFORMATION TILLBEHÖR TILL TRAKEOSTOMIKANYL



BIVONA® INNERKANYL

- Mycket flexibelt material som ger en icke-häftande yta
- Grå dragring för urtagning och montering
- Storleksmarkeringar
- Tillverkad i e-PFTE
- Rengörs med Bivona® rengöringsswab
- Innerkanyl till Bivona® modeller 6701 (TTS kuff), 60A1 (utan kuff) och 7501 (Aire-cuf®)
- Kan användas fler gånger om rengöringsinstruktionerna efterföljs
- Finns i 5 olika storlekar



Benämning	Art nr	Enhet	Avdfp
Bivona® innerkanyl, stl 7,0 mm	BRCA70	st	10
Bivona® innerkanyl, stl 7,5 mm	BRCA75	st	10
Bivona® innerkanyl, stl 8,0 mm	BRCA80	st	10
Bivona® innerkanyl, stl 8,5 mm	BRCA85	st	10
Bivona® innerkanyl, stl 9,0 mm	BRCA90	st	10
Bivona® innerkanyl, stl 7,0 mm	BRC270	st	2
Bivona® innerkanyl, stl 7,5 mm	BRC275	st	2
Bivona® innerkanyl, stl 8,0 mm	BRC280	st	2
Bivona® innerkanyl, stl 8,5 mm	BRC285	st	2
Bivona® innerkanyl, stl 9,0 mm	BRC290	st	2
Foam Cleaning Swabs	BICSWB	st	30

Smiths Medical Sverige AB

Box 1143, 164 22 Kista
Tlf. 08 - 594 772 50
Fax.08 - 594 772 69
info.sweden@smiths-medical.com
www.smiths-medical.com

Smiths Medical is part of the global technology business Smiths Group plc. Product(s) described may not be licensed or available for sale in all countries. The Smiths Medical design mark is a trademark of Smiths Medical. The symbol® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries.

©2012 Smiths Medical. All rights reserved. SE_13hGPI_BIC_1