

ETHICON

a Johnson & Johnson company



MONOCRYL*
(Poliglecaprone 25) Suture

MONOCRYL*

(Poliglecaprone 25) Monofilament

SUTURE

Description: **Monofilament** For smooth passage through tissue.
Synthetic For minimal tissue reaction.
Absorbable By hydrolysis. Essentially complete in 90 -120 days.

Tensile strength: **Post implantation** **Approximate % original strength remaining**

7 days	60% (violet) 50% (undyed)
14 days	30% (violet) 20% (undyed)
21 days	0% (undyed)
28 days	0% (violet)

Color: Violet or undyed.

Range: 6-0 to 1 (USP). Supplied as needled sutures and ligatures.

Indications: MONOCRYL* sutures are intended for use in general soft tissue approximation and/or ligation where an absorbable material is indicated. Typical areas of use include:










- Subcuticular sutures
- Small intestine anastomoses
- Urological anastomoses

Full details in the Instruction for Use included in every package.






MONOCRYL* (Poliglecaprone 25) Monofilament Suture

1/2 Circle Taper Point Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 TF	13 mm	6-0	VIOLET	70		36	Y432H
		5-0	VIOLET	70		36	Y925H
 RB-1	17 mm	5-0	VIOLET	70		36	Y213H
		4-0	VIOLET	70		36	Y214H
		4-0	VIOLET	70	VISI-BLACK*	36	Y3040H
		3-0	UNDYED	70		36	T4215AH
		3-0	VIOLET	70		36	Y215H
 RB-1 RB-1	17 mm	5-0	VIOLET	90	VISI-BLACK*	36	C3060
		4-0	VIOLET	70	VISI-BLACK*	36	Y3214H
 SH-1	22 mm	4-0	VIOLET	70		36	Y218H
		4-0	VIOLET	70	VISI-BLACK*	36	Y3100H
		3-0	VIOLET	70		36	Y219H
		3-0	VIOLET	70	VISI-BLACK*	36	Y3110H
		2-0	VIOLET	70	VISI-BLACK*	36	Y3120H
 SH-1 SH-1	22 mm	4-0	VIOLET	70	VISI-BLACK*	36	Y3218H
		4-0	VIOLET	90	VISI-BLACK*	36	Y5218H
		3-0	VIOLET	70	VISI-BLACK*	36	Y3219H
 SH	26 mm	4-0	VIOLET	70		36	Y415H
		4-0	VIOLET	70	VISI-BLACK*	36	Y3150H
		3-0	UNDYED	70		12	NP0001
		3-0	VIOLET	70		36	Y416H
		3-0	VIOLET	70	VISI-BLACK*	36	Y3160H
		3-0	VIOLET	8x45 C/R	VISI-BLACK*	12	Y3864G
		3-0	VIOLET	8X70 C/R	VISI-BLACK*	12	Y7840G
		2-0	VIOLET	70		36	Y417H
		0	VIOLET	70		36	Y418H
		2-0	VIOLET	70	VISI-BLACK*	36	Y3170H
 SH SH	26 mm	3-0	VIOLET	70	VISI-BLACK*	36	Y3416H
		3-0	VIOLET	90	VISI-BLACK*	36	Y5270H
		2-0	VIOLET	70	VISI-BLACK*	36	Y3417H
 MH-1	31 mm	3-0	VIOLET	70		12	W3200
		2-0	VIOLET	70		12	W3441
 MH	36 mm	2-0	VIOLET	70		36	Y227H

MONOCRYL* (Poliglecaprone 25) Monofilament Suture




1/2 Circle Taper Point Needles (Heavier wire diameter)

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 CT-2	26 mm	2-0	VIOLET	90		36	Y326H
		0	VIOLET	70		36	Y270H
 CT-1	36 mm	2-0	VIOLET	70		36	Y259H
		0	VIOLET	70		36	Y260H
		0	VIOLET	90		36	Y946H
		1	VIOLET	70		36	Y261H
 CTX	48 mm	2-0	VIOLET	90		12	Y979G
		1	VIOLET	90		12	Y977G

1/2 Circle ETHIGUARD* Blunt Point Needles




Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
XLHD	70 mm	1	VIOLET	90		12	W3709

5/8 Circle Taper Point Needles




Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 UR-6	26 mm	3-0	VIOLET	70		12	W3625
		2-0	VIOLET	70		36	Y605H
 UR-6	26 mm	3-0	VIOLET	70		36	Y602H
 UR-5	36 mm	2-0	VIOLET	70		36	Y375H

MONOCRYL* (Poliglecaprone 25) Monofilament Suture





1/2 Circle TAPERCUT* Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 DA-1	13 mm	5-0	UNDYED	45	Compound Curve	12	Y4950G
 V-7	26 mm	2-0	VIOLET	70		36	Y999H
 V-34	36 mm	2-0	VIOLET	70		12	W3488
		0	VIOLET	70		12	W3489
		1	VIOLET	90		36	Y919H

1/2 Circle Reverse Cutting Needles




Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 X-1	22 mm	2-0	UNDYED	70		36	C459
 CP-2	26 mm	2-0	VIOLET	70		12	W3230
 CP-1	36 mm	0	UNDYED	70		36	C267

3/8 Circle Reverse Cutting PRIME* Needles




Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 P-1	11 mm	6-0	UNDYED	45	MULTIPASS*	36	MPY489H
		6-0	UNDYED	45		36	Y489H
		5-0	UNDYED	45	MULTIPASS*	36	MPY490H
 P-3	13 mm	6-0	UNDYED	45	MULTIPASS*	36	MPY492H
		6-0	UNDYED	45		36	Y492H
		5-0	UNDYED	45	MULTIPASS*	36	MPY493H
		5-0	VIOLET	45		36	Y491H
		5-0	UNDYED	45		36	Y493H
		4-0	UNDYED	45	MULTIPASS*	36	MPY494H
		4-0	UNDYED	45		36	Y494H
 PS-3	16 mm	5-0	UNDYED	45	MULTIPASS*	36	MPY500H
		5-0	UNDYED	45		12	W3204
		4-0	UNDYED	45	MULTIPASS*	36	MPY501H
		4-0	UNDYED	45		12	W3205
 PS-4C	16 mm	5-0	UNDYED	45	Compound Curve	12	NP0127
		4-0	UNDYED	45	Compound Curve	12	NP0126

MONOCRYL* (Poliglecaprone 25) Monofilament Suture



3/8 Circle Reverse Cutting PRIME* Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 PS-2	19 mm	5-0	UNDYED	45	MULTIPASS*	36	MPY495H
		4-0	UNDYED	45	MULTIPASS*	36	MPY496H
		4-0	UNDYED	45		12	W3206
		4-0	UNDYED	70	MULTIPASS*	36	MPY426H
		3-0	UNDYED	45	MULTIPASS*	36	MPY497H
		3-0	UNDYED	45		12	W3207
		3-0	UNDYED	70		36	C427
		3-0	UNDYED	70	MULTIPASS*	36	MPY427H
 PS-1	24 mm	4-0	UNDYED	45		36	T4682H
		4-0	UNDYED	70	MULTIPASS*	36	MPY935H
		4-0	UNDYED	70		36	Y495H
		3-0	UNDYED	70		12	W3208
		3-0	UNDYED	70	MULTIPASS*	36	MPY683H
		2-0	UNDYED	70	MULTIPASS*	36	MPY684H
		2-0	UNDYED	70		36	Y684H
 PS	26 mm	4-0	UNDYED	70	MULTIPASS*	36	MPY3212H
		4-0	UNDYED	70		12	W3212
		3-0	UNDYED	70	MULTIPASS*	36	MPY3213H
		3-0	UNDYED	70		12	W3213


3/8 Circle Reverse Cutting Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 FS-2	19 mm	4-0	UNDYED	45		36	Y292H
		4-0	UNDYED	70		36	Y422H
		3-0	UNDYED	45		36	Y293H
		3-0	UNDYED	70		36	Y423H
 FS-1	24 mm	4-0	UNDYED	70		36	Y291H
		3-0	UNDYED	70		36	Y442H
		2-0	UNDYED	70		36	Y443H
 FS	26 mm	3-0	UNDYED	70		12	W3326

Straight Cutting Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 TS	48 mm	4-0	UNDYED	70		36	C1041
 KS	60 mm	4-0	UNDYED	70		12	W3649
		3-0	UNDYED	70		12	W3650

1/4 Circle Spatula Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
	8 mm	6-0	UNDYED	45	12	W3552

Ligatures

	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
LIGAPAK*	3-0	VIOLET	250	Reel	10	W3020
SUTUPAK*	2-0	VIOLET	150		12	W3113
	1	VIOLET	150		12	W3115
	0	VIOLET	150		12	W3114
SUTUPAK*	2-0	VIOLET	2x70		36	Y625H
SUTUPAK*	3-0	UNDYED	3x45		36	Y634H