



STERIKING® Daily Control G2 Test Pack

The STERIKING® second generation Daily Control test packs are single-service Bowie & Dick type tests for pre-vacuum sterilizers for steam penetration, steam quality and air leak detection. The STERIKING® Daily Control Test packs are designed to monitor the attainment of the critical variables in the Bowie & Dick test program.

Conformity to International Standards

The STERIKING® Daily Control G2 test packs conform to the following product standards: ISO 11140-1:2014 (Type 2 and 6) and ISO 11140-4:2007.

STERIKING® monitoring products are designed, validated, and manufactured to suit their intended purposes.

Technical Data & Performance Characteristics

The STERIKING® Daily Control G2 test packs are pre-assembled from approved filter materials with a non-toxic and laminated indicator sheet. The product is also latex-free.

The STERIKING® Daily Control Test packs are for daily testing of pre-vacuum steam autoclaves to demonstrate compliance with EN 285 standard and as a routine test of performance as in ISO 17665-1:2006. The indicator show a uniform dark color change after exposure to dry saturated steam at 134°C for 3,5 min. according to 11140-4:2007

Storage Recommendations & Shelf Life

It is recommended that the STERIKING® products are kept in the original, closed transport package and is stored in dry and clean conditions protected from direct sunlight and excessive moisture.

If stored according to the instructions the shelf life of STERIKING® test packs is at least three years from the date of manufacture that is marked on each test pack.

Sales and Transport Packing

STERIKING® Daily Control G2 test packs are packed in recyclable bleached corrugated cases with inner wrap. The

sales and packing unit is one case (=20 test packs). Each case bears the information required under ISO 11140-1:2014.

Please refer to the local/national regulations regarding waste disposal.

In Case of Complaint

In event of any complaint, the lot number and identification code must be provided by the complainant. For evaluation of claimed product, a defective sample (or a digital photo) and description of the defect together with an unused specimen must be made available to Wipak.

STERIKING® is a registered trademark of Wipak

Steriking® DCG2 – Daily Control Test Pack

Code	Packs/Case
DCG2	20

This specification refers to the named product group and shall be valid until the next revision. Other product related documents may be available upon request.

The information contained here is to our knowledge accurate and reliable as of the date of the publication. Wipak extends no warranties and makes no representations as to the accuracy or completeness of the information contained herein, and assumes no responsibility regarding the consequences of its use or for any printing errors. It is the customer's responsibility to inspect and test our products in order to satisfy himself as to the suitability of the products for the customer's particular purpose and suitability to the actual circumstances the product is exposed to. The customer is also responsible for the appropriate, safe and legal use, processing and handling of our products especially when recommendations for safe use and storage are given. Nothing herein shall constitute any warranty, nor is a protection from any law or patent to be inferred.