

4921**ORD[®]**
Xray
OPERATION ROOM DRESSINGS

Operation dressings

Pensos cirúrgicos

Pansements operatoires

Op-Abdecktuch

Operatieverbanden

Medicazioni per interventi

Apósitos quirúrgicos

Operations textilier

Sterile non-woven balls

- 70% viscose and 30% polyester.
- Approx. weight: 30grs/sq.m.
- X-Ray detectable thread with approx. 60% barium sulphat contents.
- The product is presented inside the cavities of the package. Each package has 5 cavities. Products are controlled by visual conference of the package and it's fulfilled five cavities.
- 5 cavities package with double control tag + Peel-open pack + shelf box + transport carton.

[click here for labeling information](#)**DIAMETER:**
DIÂMETRO:Flat size.
Diâmetro achatado.Round size.
Diâmetro arredondado.

Sterile - Double pack - Control tag



Code	REF	Weight	X Diamater	Y Diamenter	X-Ray	Gauze diameter	Pcs./Pack	Pcs./Sh.box	Pcs./Tr. carton
4921-021	1	0,27g	25	30	X	100cm2	5	250	2250
4921-022	3	0,53g	30	35	X	195cm2	5	150	1350
4921-023	4	0,92g	40	45	X	340cm2	5	150	1350
4921-024	5	1,63g	55	65	X	510cm2	5	125	750