



### Product description

Order number	77.3898.011
Product description	Multi-Safe medi 11, round disposal container 11 l, without bracket, height: 220 mm, lid?: 298 mm
UN identification	no

### Product characteristics

Cap	flat
Shape of opening	round
Removal aid - turning aid	Removal aid

### Size

Max. volume	11 l
Diameter of base	257 mm
Diameter of opening	88 mm
Height of product	220 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	yellow
Colour of cap	red
Colour of lid	yellow

### **Purity & certification**

Satisfies the requirement	TRBA 250, DIN EN ISO 23907
Product in accordance with MDD or IVDD	No marking necessary
CE certified	No CE marking
Lot or serial number	non batched

### **Packaging**

Minimum order qty.	250
Pcs. / inner box	250
Pcs. / box	250
Pcs. / pallet	250
Case volume	1,68 cbm
Weight of case	137,97 kg
EAN of inner box	4038917409483
EAN of outer box	4038917409476

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

