

4231

Finger lift strap
Tira de levantamento

White colour
Cor branca

Conf. cohes. elastic bandage

Ligadura elástica coesiva

Bande extensible cohésive

Kohäsive elastische Fixierbinde

Elastisch cohesief verband

Benda elástica aderente

Venda elástica cohesiva

Cohesive binda



Conforming cohesive elastic bandages

- 36% Cotton, 36% viscose and 28% polyamide.
- The cohesiveness is obtained by latex free glue. Non-slipping, conforming to body contours. Recommended to joints, fingers or extremities treatments, light traumatism, distortions or contusions. Also indicated for injuries caused by sport activities.
- With initial guide strap.
- Stretched Length: 4m or 20m.
- Individual shelf box + transport carton.
- Elasticity: 115%

Non sterile

Code	Size	Pcs./Sh.box	Pcs./Tr. carton
4231-001	4m x 4cm	1	72
4231-002	4m x 6cm	1	72
4231-003	4m x 8cm	1	72
4231-004	4m x 10cm	1	36
4231-020	20m x 6cm	1	16
4231-021	20m x 8cm	1	12
4231-022	20m x 10cm	1	10

