

425

Ribbon gauze

Mecha de gazeMèche de gazeTamponadestreifenGaaswiekStriscia di garzaTira de gasaTamponad av gasvöv

Gasvövstamponad

- Bomull. • För tamponering. • Uppskattad vikt: 33g/m².
- 22 trådar/cm² enligt EN 14079.
- Förpackning:
 - Osteril: Avdfp + transportkartong.
 - Steril: Peel-open pack + avdfp + transportkartong.

Osteril

Artnr	Storlek	Avdfp	Trpfp	EAN Avdp	EAN Trpfp
425-001	5m x 1,25cm	24	1728	560562240341	560562200341
425-002	5m x 2,5cm	12	864	560562240342	560562200342
425-003	5m x 5cm	12	432	560562240343	560562200343
425-004	5m x 7,5cm	8	288	560562240344	560562200344
425-005	5m x 10cm	6	216	560562240345	560562200345

Steril

Artnr	Storlek	Avdfp	Trpfp	EAN Avdp	EAN Trpfp
425-501	5m x 1,25cm	60	1080	560562241422	560562201422
425-502	5m x 2,5cm	30	540	560562241423	560562201423
425-503	5m x 5cm	26	468	560562241424	560562201424