

## ArjoWrap 366 NW



Choice of SMS non-woven for the safest packaging solution.

### Main characteristics

- Cellulose for the highest possible bacterial barrier performance
- Synthetic fibers for increased mechanical resistance

### Sterilization method

- Steam, EO

### Application

- Sterilization wrap for all tray sizes

Art. Number	Description	Specification
0153050	ArjoWrap 366 NW	50x50cm, green, 66g
0153060	ArjoWrap 366 NW	60x60cm, green, 66g
0153075	ArjoWrap 366 NW	75x75cm, green, 66g
0153090	ArjoWrap 366 NW	90x90cm, green, 66g
0153100	ArjoWrap 366 NW	100x100cm, green, 66g
0153120	ArjoWrap 366 NW	120x120cm, green, 66g