

490

Dressing set

Set de penso

Set de soins

Verband set

Verbandset

Pacchetto di medicazione

Set de cura

Omlägningsset

STERILE DRESSING SETS

All the necessary components for routine dressing procedures are grouped, placed in the order of use and sterilised in a unique pack.

The set's wrapped in a large field that, when opened, will automatically cover the working place creating immediately a sterile field. It's exclusively made by single use components.

Packing: Individual peel-open pack + shelf box + transport carton

Contents	Uni.
Wrapping field, 50x50cm	1
Tray, 18x9,5x2,5cm	1
Gauze swabs, 4x5cm	4
Gauze swabs, 7,5x7,5cm	2
Gauze balls, D	2
Green forceps	2

Code	Pcs./Pack	Pcs./Sh.box	Pcs./Tr. carton
490-002	1	50	300

Code	Pcs./Pack	Pcs./Sh.box	Pcs./Tr. carton
490-002	EAN: 560562221303	EAN: 560562241303	EAN: 560562201303