

Portex® Tracheostomy Tubes

Cuffed Blue Line Ultra® Suctionaid

Not available in USA

[INQUIRE ABOUT THIS PRODUCT](#)

PORTEX

A thermosensitive PVC tracheostomy tube with clear, soft flange and Soft Seal cuff inner cannulae available.

Incorporating "Suctionaid".

Product Features

- ▶ Blue Line Ultra® Suctionaid features an integral suction lumen to aid removal of secretions from above the cuff
- ▶ Thermosensitive PVC - provides sufficient rigidity for initial insertion, and then softens at body temperature to accommodate individual patient anatomy
- ▶ 105 degree angle for comfort in-situ
- ▶ Tube is suitably radiopaque to enable confirmation of tube position
- ▶ Soft Seal cuff - velvet soft; low pressure, high volume cuff, with larger cuff resting diameter
- ▶ Clear markings on pilot balloon provide relevant information
- ▶ Flange is soft for maximum patient comfort, and clear to ensure aesthetic acceptability
- ▶ Obturator features special clip design to minimise tube tip movement during insertion
- ▶ Supplied with tapes



Ordering Information:

Box quantity: 10



[Symbol Key](#)

[Product Codes](#)

[Documents](#)

[Associated Products](#)

Product Specifications



Cuffed Blue Line Ultra® "Suctionaid" Tracheostomy Tube

Reorder Code	ID (mm)
100/860/060	6.0mm
100/860/070	7.0mm
100/860/075	7.5mm
100/860/080	8.0mm
100/860/085	8.5mm
100/860/090	9.0mm
100/860/100	10.0mm

Notes: Inner cannulae available - supplied in packs of 2 units. Plain inner cannulae ref: 100/850