

OXYGEN THERAPY

Oxygen Masks



13113

Description

- Soft anatomical form ensures good fit and patient comfort
- Adjustable metal nose piece for security and comfort
- Wide adjustable elastic strap with sealed ends fits securely and comfortably below the ears

REF	DESCRIPTION	QTY
13110^a	Salter-Style™ Adult oxygen mask with 213cm crush resistant oxygen supply tube, elastic strap style	50
13113^a	Salter-Style™ Paediatric oxygen mask with 213cm crush resistant oxygen supply tube, elastic strap style	50
221/5001^a	DAR™ Adult oxygen mask with 200cm crush resistant oxygen supply tube, elastic strap style	25
211/5013^a	DAR™ Paediatric oxygen mask with 200cm crush resistant oxygen supply tube, elastic strap style	25

^a Not made with DEHP