


# Surshield™ Versatus™ -WP

Winged and Ported I.V. catheter with Protector



## Product Information



Description	Raw Materials
1 Catheter	Polyurethane + BaSO <sub>4</sub>
2 Winged Catheter hub Caulking pin	Polypropylene Polyoxymethylene
3 Cannula (internal needle)	Stainless steel
4 Needle hub	Acrylonitrile butadiene styrene
5 Luer-lock cap	Polypropylene
6 Vent fitting filter Vent fitting body	Polyethylene Acrylonitrile butadiene styrene
7 Port Closure	Polypropylene + master batch
8 Injection port valve	Silicone tube
Protector	Polyethylene
<b>Surshield Protector</b>	
9 Protection housing	Polyamide
Safety fingers	Polysulfone
Sliding bush + Needle blocker	Polycarbonate + stainless steel cup

Does not contain components made of PVC or natural rubber latex.

### Developed and manufactured by skilled technicians in Europe

The Surshield protector allows to completely close the needle, in order to prevent accidental needle stick injuries when disposing of the intravenous cannula. Using a safety catheter does not require any change in the insertion technique.

#### Sterilization

With ethylene oxide to a Sterility Assurance Level of at least 10<sup>-6</sup>

#### Storage

- Shelf life of 5 years
- The batch number and expiry date are printed on the paper of the blister pack
- Do not store at extreme temperature and humidity

#### Precautions

- Do not use if unit package is damaged
- Do not re-use / do not re-sterilize
- Do not attempt to re-insert a partially or completely withdrawn needle
- Dispose of safely after single use to avoid risk of infection
- Application time should be consistent with accepted standards of practice

# Surshield™ Versatus™ -WP

Winged and Ported I.V. catheter with Protector



## Your catheters of choice, part of a wider IV portfolio.

### Product Range

Product Code	Colour	External Catheter Ø	Internal Catheter Ø	Catheter Length	Flow Rate ml/min.
SR+DS2619PD	Violet	26G 0.64 mm	0.45 mm	19 mm	19
SR+DS2419PD	Yellow	24G 0.74 mm	0.55 mm	19 mm	29
SR+DS2225PD	Blue	22G 0.90 mm	0.65 mm	25 mm	42
SR+DS2032PD	Pink	20G 1.00 mm	0.75 mm	32 mm	59
SR+DS1832PD	Green	18G 1.30 mm	0.95 mm	32 mm	103
SR+DS1845PD	Green	18G 1.30 mm	0.95 mm	45 mm	96
SR+DS1745PD	White	17G 1.50 mm	1.15 mm	45 mm	155
SR+DS1645PD	Grey	16G 1.75 mm	1.35 mm	45 mm	225
SR+DS1445PD	Orange	14G 2.00 mm	1.55 mm	45 mm	290



### Packaging

**Single Unit** Blister pack

**Unit Box** 50 pieces

**Shipping Carton** 200 pieces  
Dimensions 390 x 384 x 180 mm

### Regulatory Status

CE Certification – Class IIa medical device

#### Notified body

Istituto Superiore di Sanità 0373

#### Manufacturer

Delta Med S.p.A. (Italy)

#### Distributor

TERUMO EUROPE N.V.,  
Interleuvenlaan 40, 3001 Leuven, Belgium

This range complements the wider Terumo product offering including a comprehensive portfolio of both safety and conventional IV products which can be found at our web page on [www.terumo-europe.com](http://www.terumo-europe.com)



Terumo Corporation  
+81 3 6742 8500

Terumo Europe N.V.  
+32 16 38 12 11

Your Terumo partner