


Spirometry

Each patient needs specific optimized ventilation during anesthesia and respiratory care. Our patient spirometry provides continuous and non-invasive measurements of real patient values for pulmonary mechanics and ventilation.



733910-HEL



733950

Spirometry Kits and Accessories		
Product Number	Description	Case Quantity
8002718	Disposable patient spirometry kit, pediatric <ul style="list-style-type: none"> • One anesthesia gas sampling line, 7 ft (2 m) • One Pedi-lite+ flow sensor • One spirometry tube, 7 ft (2 m) <i>All parts are pre-connected, use with CAiOV, CAiOVX, COV, COVX or Capnomac Ultra*</i>	50
889560	Disposable patient spirometry kit, adult <ul style="list-style-type: none"> • One anesthesia gas sampling line, 10 ft (3 m) • One D-lite* sensor, yellow • One spirometry tube, yellow, 10 ft (3 m) <i>All parts are pre-connected, for use with CAiOV and Capnomac Ultra</i>	50
M1032634	Disposable patient spirometry kit, pediatric <ul style="list-style-type: none"> • One anesthesia gas sampling line, 10 ft (3 m) • One single-use Ped-lite+ sensor • One spirometry tube, 10 ft (3 m) 	5
883874	Gas and spirometry accessory kit <ul style="list-style-type: none"> • 10 disposable anesthesia gas sampling lines, 10 ft (3 m) (73319) • 10 single-use straight T-adapters (73385) • 1 D-fend water trap (876446) • 2 spirometry tubes • 2 reusable D-lite sensors (733910) 	1
M1145081†	Patient spirometry kit, adult	50
M1145077†	Disposable D-lite sensor, yellow	50
733910-HEL	Reusable D-lite sensor, 22 mm M, 15 mm F, 15 mm M; material: transparent yellow polyphenylsulfone; dead space: 9.5 ml; resistance: 0.5 cm H ₂ O at 30 l/min; patient range: TV 150 - 2000 ml	1
73393	Reusable Ped-lite+ sensor, 22 mm M, 15 mm F, 15 mm M; material: transparent yellow polyphenylsulfone; dead space: 2.5 ml; resistance: 1.0 cm H ₂ O at 10 l/min; patient range: TV 15 - 300 ml	1
733950	Disposable D-lite sensor, 22 mm M, 15 mm F, 15 mm M; material: yellow polystyrene; dead space: 9.5 ml; resistance: 0.5 cm H ₂ O at 30 l/min; patient range: TV 150 - 2000 ml	50
884101	Spirometry tube, yellow, 10 ft (3 m)	5
890031	Spirometry tube, yellow, 7 ft (2 m)	5
891191-HEL	Spirometry tube, yellow, 20 ft (6 m)	5
M1145088†	Disposable spirometry tube, yellow, 10 ft (3 m)	5
M1145087†	Disposable spirometry tube, yellow, 7 ft (2 m)	5
M1145089†	Disposable spirometry tube, yellow, 20 ft (6 m)	5

† Product is only available in the U.K.

M - male F - female TV - tidal volume