

Alaris™ GP Series Oncology primary set



Alaris™ GP Series Oncology primary set 60950E

- Needle-Free Connectors
- Set Length 260 cm
- Priming Volume 24 ml
- Bore 3 mm
- 1 Spike(s) (Vented/Unvented)
- In-Line Drip Chamber
- Drops per ml 20
- Hydrophobic Air Inlet Filter (with Cap)
- 15 µm Drip Chamber Filter
- 5 SmartSite® Valve Y Port
- 2 Pinch Clamp(s)
- 1 Roller Clamp(s)
- Alaris™ Safety Clamp
- Pump Segment
- 1 Male Luer(s)
- Rotating Male Luer(s)
- Hydrophobic Priming Cap
- Polyvinyl Chloride (PVC)
- Natural rubber latex is not part of the material formulation
- DEHP is not part of the material formulation
- Quantity 50



5 SmartSite® 1 Male Luer(s)
Valve Y Port

0086